SECTION 1 - NO POWER

Symptoms:

- The batteries are fully charged.
- All electrical components are connected correctly.
- The on/off key is pressed and the power does not come on.

Diagnosis:

The power has been interrupted somewhere in the system.

Solution:

Use the following procedure to find the source of the interruption:

- 1. Measure voltage across pin 1 (B+) and pin 2 (B-) on connector 1a. See figure 2.
- If your multimeter indicates OVDC, then go to the next step.
- *If your multimeter indicates about 25VDC*, then replace the VSI controller (1) and retest the system.
- If your multimeter indicates 0 18VDC, then recharge the batteries and retest the system.
- 2. Unplug connector 1c from connector 2a. See diagram 2.
- 3. Measure voltage across pin 6 (B+) and pin 8 (B-) on connector 2a and across pin 7 (B-) and pin 5 (B+) on connector 2a. See figure 3.
- *If your multimeter indicates less than 18VDC*, then go to the next step.
- *If your multimeter indicates more than 18VDC*, then replace the VSI controller (1) and retest the system.
- 4. Remove the seat and foot platform assembly. Refer to the power base owner's manual.
- 5. Remove the shroud. See figure 4.
- 6. Measure voltage across connector 6b and connector 7c. (*If your multimeter indicates 0VDC, then measure voltage across connectors 6c and 7b.*) See diagram 2.
- If your multimeter indicates 0VDC or more than 18VDC, then go to the next step.
- If your multimeter indicates 0VDC 18VDC, then recharge the batteries and retest the system.
- If the batteries do not appear to be taking a charge, then go to "Flash Code #1 - Low Battery", step 16.

WARNING! Battery posts, terminals, and related accessories contain lead and lead compounds. Wash hands after handling.

WARNING! Always protect the batteries from freezing and never charge a frozen battery. Charging a frozen battery may result in personal injury and/or damage to the battery.




Figure 2. Connector 1a




Figure 3. Connector 2a




Figure 4. Jazzy 1121 Shroud Assembly/Disassembly



MANDATORY! RED (+) cables must be connected to positive (+) battery terminals/ posts. BLACK (-) cables must be connected to negative (-) battery terminals/posts. Failure to connect the battery cables and harnesses in the proper manner may result in personal injury and/or damage to the power chair. REPLACE cables immediately if damaged.

- 7. Unplug connector 7a from connector 2b and connector 6a from 2c. See diagram 2.
- 8. Measure voltage across connector 6b and connector 6c. See diagram 2.
- 9. Measure voltage across pin 1 and pin 2 on connector 6a. See figure 5.
- If your multimeter indicates the same voltage as measured in step 8, then go to the next step.
- If your multimeter indicates a different voltage (by at least 0.2VDC), then replace the rear battery harness (6) and retest the system.
- 10. Measure voltage across connector 7b and connector 7c. See diagram 2.
- Measure voltage across pin 1 and pin 2 on connector 7a.
 See figure 5.
- If your multimeter indicates the same voltage as measured in step 10, then go to the next step.
- If your multimeter indicates a different voltage (by at least 0.2VDC), then replace the front battery harness (7) and retest the system.
- 12. Unscrew the electronics tray assembly and lift it off of the 1121 frame. **See figure 6.**
- 13. Measure resistance across contact 1 and contact 2 of the circuit breaker (2f). See figure 7.
- If your multimeter indicates less than 1 ohm, then replace the power interface harness (2) and retest the system.
- *If your multimeter indicates an open,* then replace the circuit breaker (2f) and retest the system.




Figure 5. Connector 6a (also connector 7a)




Figure 6. Jazzy 1121 ElectronicsTray Assembly/ Disassembly




Figure 7. Circuit Breaker (2f)