



Allied Medical

Short Term Loan and

Rehabilitation Solutions



**ALLIED  
MEDICAL**  
LIMITED

# Product Range

## MANUAL WHEELCHAIRS 04

ADULT WHEELCHAIRS 04

PAEDIATRIC WHEELCHAIRS 07

BARIATRIC WHEELCHAIRS 08

WHEELCHAIR CUSHIONS 09

WHEELCHAIR ACCESSORIES 11

## POWERCHAIRS 14

## HOME HEALTHCARE 15

BATHING 15

TOILET AIDS 17

LIVING AIDS 19

## WALKING AIDS 26

ROLLATORS 26

CANES 28

CRUTCHES 31

## BRACING & SUPPORTS 33

## EXERCISE PRODUCTS 51

## PAIN MANAGEMENT 52



## Karma Soma Agile Wheelchair

**KMAG1616P 16"x 16" Seat**

**KMAG1818P 18"x 18" Seat**

**KMAG1816P 18"x 16" Seat**

**KMAG2018P 20"x 18" Seat**

- Standard flip-back armrests
- 80° removable footrests to ensure safe transfer in and out of the wheelchair
- Compact and lightweight for easy transport
- Quick-release rear wheels for lighter carrying weight and smaller storage space
- Built with aircraft-grade aluminium
- Elevating legrests and anti-tippers available on request
- Transit kit available on request
- Pearl Silver or black
- Unit Weight 8.8 kg
- Maximum User Weight 100 kg



## Karma S-Ergo Self Propelling

**KM125Q2416 16"x 17" Seat**

**KM125Q2418 18"x 17" Seat**

- Flip-back, ergonomically shaped armrests designed to offer comfort and support
- Extended armrests allow for wider seat area and minimise the overall width
- Ergo push rim that aids wheelchair propulsion
- Swing-away footrests swing both inwards & outwards
- Better side support of legs from "tube-in-centre" footplate design
- Aegis upholstery anti-bacterial shielding
- Optional anti-tippers
- Quick release 24" rear wheels
- Weighs under 13.2 kg
- Maximum User Weight 115 kg



## Karma Lightweight Wheelchair Self Propelling or Transit

**KM8020QDXX (Self Propelling)**

**KM8520FXX (Transit)**

- Versatile and easily adjustable - available in 7 sizes
- Lightweight aircraft grade aluminium frame
- Folding frame with reinforced cross bracing
- Removable rear wheels on self propelling model
- Flip back, height adjustable armrests: desk or full length arm pad
- Height adjustable backrest
- Push to lock brakes
- Pelvic positioning belt
- 70-degree legrests with calf strap
- Multiple-position rear wheel mount
- Unit Weight 15.5 kg
- Maximum User Weight 130 kg (14" - 18" width frames), 160 kg (20" - 22" width frames)







## Karma Flexx Wheelchair

**KMFLEXXSPXX (Self Propelling)**

**KMFLEXXSTRXX (Transit)**

**Seat Width from 14" to 18" and Depth 16" to 18"**

- Fixed angle adjustment in back canes -4 to +16° in set increments
- 15.5" minimum seat to floor heights (with 22" wheels and 5" castors)
- Rise and fall in the front castors and 25 rear axle positions
- 80° swing in / swing out legrest hangers and flip back height adjustable armrests
- Options including elevating legrests, 90° hangers, stump supports and attendant brakes available
- Maximum User Weight 130kg



## Karma Flexx HD Wheelchair

**KMFLEXXHDSPXX (Self Propelling)**

**KMFLEXXHDTRXX Transit)**

**Seat Width from 20" to 24" and Depth 18" to 20"**

- Fixed angle adjustment in back canes -4 to +16° in set increments
- Rise and fall in the front castors and 25 rear axle positions
- 80° swing in / swing out legrest hangers with aluminium footplates
- Flip back height adjustable armrests
- Options including elevating legrests, 90° hangers, stump supports and attendant brakes available
- Maximum User Weight 170 kg
- Tie downs are not available



## Karma Standard Lightweight Self-Propelling Reclining Wheelchair

**KM5000-16 (16"x 16" Seat)**

**KM5000-20 (20"x 18" Seat)**

**KM5000-18 (18"x 18" Seat)**

- Reclining backrest
- Headrest with movable cushion to support the head
- Quick release, removable rear wheels for transport
- Fully removable full length armrests
- High mounted, push-to-lock brakes
- Elevating swing-away leg rests, with calf pad
- Reinforced cross bracing – strong but easily folded
- Multiple position rear wheel and castor mounts to optimise centre of gravity and seat to floor height adjustments
- Maximum User Weight 100 kg





## AML Self Propelling Wheelchair

**AMLWCSP16 (16" x 16" Seat) AMLWCSP18 (18" x 16" Seat)**

- Powder coated steel frame
- Detachable full armrest
- Swing-away footrest
- Steel chrome handle rim
- Plastic black footplate
- Plastic front fork
- Folds for easy storage and transport
- PU arm pads for extra comfort
- Nylon upholstery
- 24" real spoke wheels with PU tires
- 8" x 1.25" front castors with PU tires
- Unit Weight 19.5 kg
- Maximum User Weight 135 kg



## AML Transit Wheelchair

**AMLWCTR16 (16" x 16" Seat) AMLWCTR18 (18" x 16" Seat)**

- Powder coated steel frame
- Detachable full armrest
- Swing-away footrest
- Steel chrome handle rim
- Plastic black footplate
- Plastic front fork
- Folds for easy storage and transport
- PU arm pads for extra comfort
- Nylon upholstery
- 12.5" PU mag wheels
- 8" x 1.25" front castors with PU tires
- Unit Weight 19.5 kg
- Maximum User Weight 150 kg



## Karma Flexx Junior Self Propelling Wheelchair

**KMFLEXXSPJR1011 (10" x 11" Seat)**

**KMFLEXXSPJR1212 (12" x 12" Seat)**

**KMFLEXXSPJR1414 (14" x 14" Seat)**

- Configurable dimensions
- Height adjustable flip-back armrests
- Swing in and out 90° legrest, adjustable length
- Padded height and angle adjustable backrest
- Reinforced double crossbrace
- Quick release pneumatic rear wheels
- Push to lock brakes with brake extension
- 2 point lap belt
- Heel loops
- Maximum User Weight limit of up to 75 kg



## Karma Reclining Paediatric Self Propelling Wheelchair

**KM7520R13**

- Easily adjustable push handle to accommodate carers of all heights
- Quick release, removable rear wheels
- Fully removable, height adjustable armrests
- High mounted, push to lock brakes, that can be relocated in numerous positions
- Swingaway 90 degree legrests with calf strap
- Reinforced cross bracing - strong but easily folded
- Durable upholstery, with hidden screws
- Multiple-position rear wheel and castor mounts to optimise centre of gravity and seat-to-floor height adjustments
- Pump and tool kit for ongoing standard maintenance



## Karma Lightweight Paediatric Self Propelling Wheelchair

**KM752011 (11.5" x 11.5" Seat)**

**KM752013 (13.5" x 13.5" Seat)**

- Easily adjustable push handle
- Quick release, removable rear wheels for ease of storing in vehicle
- Fully removable, height adjustable armrests, that allow for easier transfers
- High mounted, push to lock brakes, that can be set to numerous positions
- Swingaway 90° legrests with calf strap
- Reinforced cross bracing - strong but easily folded
- Durable upholstery with hidden screws
- Multiple-position rear wheel and castor mounts Pump and tool kit for ongoing standard maintenance



## Karma BT10 Bariatric Wheelchair

**KMBT102220 (22" x 20" Seat)**

**KMBT102420 (24" x 20" Seat)**

**KMBT102620 (26" x 20" Seat)**

- Angle adjustable backrest (from -5 degrees to +15 degrees)
- Fully detachable flip-back armrests
- Press down, swing-away detachable leg rests with angle adjustable footplates
- Rigidising bar to add strength to frame
- Swing-away anti-tippers
- Hip belt
- Calf strap
- Dual brake system – attendant brake as well as manual brakes
- Pump and tool kit for ongoing standard maintenance
- Maximum User Weight 250 kg



## AML Self Propelling Bariatric Wheelchair

**AMLWCSPB20 (20" x 18" Seat)**

**AMLWCSPB22 (22" x 18" Seat)**

**AMLWCSPB24 (24" x 18" Seat)**

- Chrome polish finish steel frame
- Steel black side-guards
- Detachable full armrest
- Padded armpads for extra comfort
- Aluminium swing-away footrest
- Folds for easy storage
- Padded PVC upholstery
- 24" mag wheels with PU tires
- 8" x 2" front castors with PU tires
- Maximum User Weight 250 kg



## AML Transit Bariatric Wheelchair

**AMLWCTRB22**

- Powder-coated steel frame with double crossbar
- Detachable full armrest with PU armpads for extra comfort
- Swing-away footrest
- Plastic black footplate
- Aluminium front fork
- Folds for easy storage and transport
- Padded nylon upholstery
- 12" rear wheels with PU tyres
- 8" x 2" front castors with PU tyres
- Maximum User Weight 230 kg





## Varilite Reflex Wheelchair Cushion

- VL73500114162 14"x 16" Seat
- VL73500116162 16"x 16" Seat
- VL73500116182 16"x 18" Seat
- VL73500118162 18"x 16" Seat
- VL73500118182 18"x 18" Seat
- VL73500118202 18"x 20" Seat
- VL73500120182 20"x 18" Seat
- VL73500120202 20"x 20" Seat

- Non-adjustable skin protection cushion featuring VARILITE's unique Air-Release Device
- Achieves adjustment through a self-regulating air chamber within the cushion that adjusts the contour of the cushion to the individual sitting on it
- Ultra-lightweight construction
- Removable incontinent cover that can be machine-washed when the need arises
- Maximum User Weight 295 kg



## Varilite Stratus Wheelchair Cushion

- VL56640 (16X16 INCONTINENT)
- VL56840 (16X18 INCONTINENT)
- VL58640 (18X16 INCONTINENT)
- VL58840 (18X18 INCONTINENT)

- For moderate to low risk of tissue breakdown and moderate positioning needs
- Unique self-inflating valve, which allows air in and out of the cushion so it can contour to the user's shape
- Two types of foam create the functional areas of the cushion: medium foam for the thigh troughs and pelvic bucketing; and firm foam for the perimeter and the medial thigh separator
- User's weight is distributed evenly over the surface of the cushion
- Maximum User Weight 320 kg



## Varilite Evolution Wheelchair Cushion

- VL74620 (14" x 16" Seat)
- VL74820 (14" x 18" Seat)
- VL76620 (16" x 16" Seat)
- VL76820 (16" x 18" Seat)
- VL76020 (16" x 20" Seat)
- VL78620 (18" x 16" Seat)
- VL78820 (18" x 18" Seat)
- VL78020 (18" x 20" Seat)

- Constructed for users at high risk of tissue breakdown and who have moderate positioning needs
- For motor and neurological dysfunction due to spinal cord injury, stroke, multiple sclerosis, cerebral palsy, and traumatic brain injury are the most frequent, though not the only, users of the Evolution cushion
- We recommend purchasing a spare cover with all cushion purchases



## Ki Mobility Axiom G Cushion

- |                          |                          |
|--------------------------|--------------------------|
| KIXG1414 (14"x 14" Seat) | KIXG1816 (18"x 16" Seat) |
| KIXG1616 (16"x 16" Seat) | KIXG1818 (18"x 18" Seat) |
| KIXG1618 (16"x 18" Seat) | KIXG2018 (20"x 18" Seat) |

- Mild contoured well supports the trochanters, transitioning the femurs and enveloping the ischial tuberosities and coccyx for support and stability
- Pre-Contoured Shape and Enhanced Positioning
- Surface Tension Reduction
- High-Quality Material Feature Rich Cover



## AML Wheelchair Seat Cushion

- |                            |                            |
|----------------------------|----------------------------|
| AMLWSC1616 (16"x 16" Seat) | AMLWSC1818 (18"x 18" Seat) |
| AMLWSC1816 (18"x 16" Seat) | AMLWSC2018 (20"x 18" Seat) |

- Ergonomically designed PU foam contoured cushion for wheelchair seats



## AML Wheelchair Anti-tippers - Pair

### AMLWCAT

- Replacement pair of anti-tippers for the AML standard self-propelling and transit wheelchairs
- Sold as a pair



## AML Wheelchair Bariatric Anti-tippers

### AMLWCATB

- Replacement pair of Anti-Tippers for the AML Bariatric wheelchairs
- Sold as a pair



## AML Wheelchair Elevating Legrest

### AMLWCELRL (LHS)

### AMLWCELRR (RHS)

- Padded Elevating Legrest for the AML wheelchairs
- Ideal for short term use while leg is in plaster or needs elevation
- Sold individually (left or right)



## Ottobock Channel Forearm Pads

### OB436A2=3-7 (Small)

### OB436A2=2-7 (Medium)

### OB436A2=1-7 (Large)

- Arm gutter to help position arm on wheelchair armrest
- Provides functional support for upper extremity weakness or paralysis
- Mounted directly onto the armrest hardware
- Fits either left or right side





## Ottobock Channel Forearm Pad

**OB436A3=2-7**

- One piece channel forearm pad
- Provide functional support if you have upper extremity weakness or paralysis
- It can be mounted directly onto the armrest hardware and can also be used with a Swivel Unit or Offset Plate to get the correct position



## Ottobock Modular Hand Pad - Horn Type

**OB436P3=R2-7 (Medium Right)**

**OB436P3=R1-7 (Large Right)**

**OB436P3=L2-7 (Medium Left)**

**OB436P3=L1-7 (Large Left)**

- Can be mounted in connection with the Ottobock modular Channel Forearm Pads
- For people with hand contractures as it can assist in preventing the hand from completely clenching



## Ottobock Modular Hand Pad - Cone Type

**OB436P4=R2-7 (Medium Right)**

**OB436P4=R1-7 (Large Right)**

**OB436P4=L2-7 (Medium Left)**

**OB436P4=L1-7 (Large Left)**

- Can be mounted in connection with the Ottobock modular Channel Forearm Pads
- For people with hand contractures as it can assist in preventing the hand from completely clenching
- Constructed of polyurethane foam



## Ottobock Modular Hand Pad - Extensor Concave

**OB436P2=R2-7 (Medium Right)**

**OB436P2=R1-7 (Large Right)**

**OB436P2=L2-7 (Medium Left)**

**OB436P2=L1-7 (Large Left)**

- Ideal for keeping the hand in a relatively neutral position
- Support surface has a built-up palmar arch for mild to moderate positioning needs
- Constructed of polyurethane foam.



## Ottobock Modular Hand Pad - Flat

**OB436P5=2-7 (Medium)**

**OB436P5=1-7 (Large)**

- If no special type of hand pad is required, the flat pad may be a possible solution
- These can also be used for the basis of a custom fabrication
- Constructed of polyurethane foam



## Gel Ovation Arm Pad

**GO310G (Size 3.5"x10")**

**GO316G (Size 3.5"x16")**

**GO312G (Size 3.5"x12")**

**GO212G (Size 2"x12")**

**GO314G (Size 3.5"x14")**

**GO214G (Size 2"x14")**

- Easily bolts onto existing armrest T sections
- They consist of a solid wood base gel pad with a neoprene nylon cover and plastic base guard
- Synergel Gel does not absorb liquids and remains consistently soft regardless of ambient temperature.



## Karma KP25.2 Powerchair

**KP2521616L 16" x 16" Seat**

**KP2521818L 18"x 18" Seat**

**KP2522018L 20"x 18" Seat**

- The new Karma KP25.2 is the updated version of the very durable KP25
- It is a cost effective rear-wheel drive powerchair that provides the user with independence and mobility
- The frame is lightweight, compact and easy to assemble. It can also be folded down to fit into a car boot
- The KP25.2 features a VR2 controller



## Jazzy 623 2.0

**PRIDEJ6232PB**

- Compact full-size mid-wheel drive powerchair with six wheels on the ground for maximum stability
- VR2 controller allowing up to two powered functions from Power tilt, Power elevate and Power elevating leg rests, swing-away or centre mounted
- High efficiency and mileage with Linux 2 pole motor
- Compatible with TRU-Balance Power Positioning
- Quick Height Adjustable (QHA) armrests
- Lower seat to floor heights (down to 17.5" with tilt)
- Maximum User Weight 136 kg
- Attendant control





## Shower Stool With Arms

### MP20106

- Lightweight and easy to move it can be used in any room
- The powder coated steel frame has four adjustable height legs and also allows for an adjustable seat angle
- Maximum User Weight 120 kg



## Height Adjustable Shower Stool with Arms

### AML20107

- Height adjustable shower stool with injection moulded plastic seat
- Powder coated frame
- Armrests assist safe transfer



## Height Adjustable Shower Chair With Arms

### AML20108

- Height adjustable shower chair with injection moulded plastic seat
- Features a curved plastic seat for comfort and drainage holes for excess water
- Armrests assist safe transfer



## Bath Board

### MP20301 (Medium)

### AML20303 (Large)

- Designed to provide the user with a comfortable and sturdy seat for independent bathing
- One piece injection moulded bath board
- Has four adjustment knobs which enable the bath board to be fitted into the curvature of the bath
- Maximum User Weight 250 kg



## Plastic Bath Seat with Suction Cups

### MP20302

- This bath seat provides a stable sitting platform that fits snugly on the bottom of the bath
- It features a cut away area which allows for easier personal cleaning
- It fixes in place for four strong suction cups
- Maximum User Weight 190 kg





## 4" Raised Toilet Seat With Armrests

### MP20402

- A moulded plastic 4" Raised Toilet Seat with padded arms that allow for security and independence when getting on and off the toilet
- Ideal for use after hip or knee replacement surgery
- It easily fits to the toilet bowl via adjustable brackets
- Maximum User Weight 190 kg



## 4" Raised Toilet Seat Without Armrests

### MP20404

- A one-piece injection moulded 4" high raised toilet seat that fits easily to the toilet bowl with three adjustable brackets
- Maximum User Weight 190 kg



## 6" Raised Toilet Seat

### DM12066

- Designed to be used on an existing toilet, it easily attaches to the toilet bowl and locks into place with two rear locks
- There are two hygiene cut-outs for ease of use and with no lid attached to the seat, it is very easy to clean
- No tools are required for installation
- Maximum User Weight 200 kg



## Bedside Commode

### MP20502

- Padded seat & back, moulded armrests
- Powder coated blue steel frame
- Height adjustable legs
- Maximum User Weight 125 kg



## 3 in 1 Commode

### AML20505

- Multifunctional - can be used as bedside commode, toilet safety frame and raised toilet seat.
- Removable backrest
- Comfortable, stable and lightweight
- Powder coated steel frame
- Comes complete with bucket, lid and splash guard



## Belgrave Bariatric Commode

### DM190

- A robust height adjustable commode with a non-peel frame and a cellular moulded backrest and over seat that can be removed when needed
- Bucket holds 14 litres and comes with a lid
- Arms can also be removed for easy side lateral transfer
- Maximum User Weight 280 kg



## Toilet Safety Frame

### MP20703

- Powder Coated Steel frame with plastic arm rests that fits around all three sides of a toilet
- Ideal for use during recovery from hip or knee surgery to ensure extra stability is provided
- The height adjustable legs have rubber tips that ensure it does not slide when being used
- Maximum User Weight 130 kg





## Good Grips Jar Opener

**GG39211**

- Open any size jar effortlessly as the sharp stainless steel teeth grips the lid firmly
- Simply hold the Jar Opener by the handle, wedge onto the jar lid, and twist to open
- The soft, non-slip handle cushions your hand and absorbs pressure
- Dishwasher safe



## AML Kettle Pourer

**SLHA4270**

- The frame helps prevent spillage
- Made of wire frame with a cradle
- Holds most kettles and pivots on a raised wire frame to allow gentle pressure to be applied to the handle to pour into the cup
- The kettle sits in the cradle and is secured by a strap



## AML Four Wheeled Walking Trolley

**AMLWALKT**

- Height adjustable with a 5" height adjustment range
- The walker is lightweight at only 6 kg.
- High strength steel tubing
- Two trays, capable of holding 15 kg
- Four wheels available for easy mobility with limited pushing



## Over Bed/Chair Table

**MP30203**

- Adjustable tilting table to be used by someone in bed or sitting in a chair
- Walnut grained melamine laminate for easy care
- Has a solid rim to prevent things from sliding off
- Fitted with four 50 mm smooth glide castors





## Perching Stool with Arms and Backrest

### AML30305

- The Perching Stool offers a padded seat and backrest with removable curved armrests that assist with safe transfers
- The rubber ferrules give a sturdy grip on the floor and the chair is height adjustable to suit a variety of users
- 4 options in 1 frame



## Perching Stool

### MP30304

- Designed to allow the user to rest and then stand again with ease, which works well for those who tire easily or recovering from surgery or an injury
- Rigid steel, powder coated frame that is able to be used in the kitchen or bathroom
- Maximum User Weight 125 kg



## AML Bed Cradle

### AMLBCA

- White powder coated steel frame bed cradle helping to eliminate the weight of blankets on feet and ankles
- 5 height adjustments
- Lower frame slides underneath
- Light weight



## KCare Patient Assist Rail

### KC9077

- Rail can be used with any sized bed, simply slide it under the mattress for easy application
- 45 mm gap between the shield and the bar, allowing for a safe and comfortable grip
- Enclosed return to minimise injury risk
- Inner-shield preventing entrapment risk
- Powder-coated tubular steel frame material



## Parsons Comfort Grip Cutlery Individual Cutlery

**PA16T001RD (Knife)**

**PA16T002RD (Fork)**

**PA16T003RD (Teaspoon)**

**PA16T004RD (Soup spoon)**

- Get maximum user control with minimum effort as the soft vinyl coating has finger indentations for a secure grip
- Also available as full set



## Parsons Comfort Grip Weighted Individual Cutlery

**PA16T400 (Knife)**

**PA16T401 (Fork)**

**PA16T403 (Soup spoon)**

- Weighted Cutlery for persons with Parkinson's or other impairments resulting in hand tremors
- Latex-free
- Home dishwasher safe
- Straight handled universal option only
- Cutlery each sold individually
- Right and left-handed versions are available for purchase by indent order only



## AML Utensil Strap

**SLHA4196**

- Ideal for people with a limited grip
- Simply attach the strap to the utensil and slip a hand through the loop
- It fits almost any hand. This item does not include the spoon
- This utensil strap is best suited for the AML utensils with large grips and looped ends



## AML Large Scoop Bowl

**SLHA4248**

- The large mould side to this scoop bowl accommodates individuals who have difficulty when eating particularly with one hand
- It has a non-slip suction cup to keep the bowl in place
- Microwave and household dishwasher safe up to 70° C except for the suction base
- White only



## AML Scoop Dish

**SLHA4249**

- The high back on the scoop dish accommodates individuals who have difficulty when eating particularly with one hand
- It has a non-slip suction cup to keep the dish in place
- Microwave and household dishwasher safe up to 70° C except for the suction base
- White only



## AML Plastic Plate Guard

**SLHA4250**

- An easy to fit and remove plastic plate guard that fits tightly around the rim of any circular plate between 200-300 mm in diameter
- The edge is curved to prevent food from falling off the plate, particularly if eating with one hand
- Microwave and household dishwasher safe up to 70° C
- White only



## Parsons Inner Lip Plate

**PA16T118 (Sandstone)**

**PA16T118BL (Blue)**

- Holds food on the plate while the user brings the fork or spoon to the edge of the plate making it ideally suited for people with use of only one hand.
- Break proof
- Microwavable
- Available in Blue or Sandstone



## Parsons Double Handled Mug With Lid

PA16T123

- The soft cushion grip is comfortable and easy to hold
- Flexible section to allow the fork to be bent to any angle
- Household dishwasher safe up to 90° C
- Large soft cushion grip to hold easily
- Features a loop at the base for attaching a utensil strap
- Available in Red and White



## AML Transparent Drinking Cup

SLHA4262

- Transparent cup with a drinking lid to help prevent spillage
- Two lids of varying sizes (large hole and small hole), allowing the user to choose the lid that makes drinking easier
- Microwave and household dishwasher safe up to 70°C



## AML Dressing Stick

SLDA5201


- A lightweight device for positioning and locating clothing featuring a smooth wooden shaft hook with a large plastic coated "C" hook at one end and a specially shaped "push/pull" hook at the other
- Length is 685 mm



## Sock and Stocking Aid

MP30805

- Lightweight, terry covered plastic cover sock aid makes it easy for socks to be put on, even if you can bend down
- The terry on the outside helps hold the sock in place and the slits in the plastic allow for it to bend around the heel
- The pulling strap makes it possible to use with just one hand
- Fully washable



## Parsons Standard Stocking Hose Aid

### PA16D005

- A flexible sock aid with nylon webbing straps
- The curved shape holds the stocking on the aid while pulling it up your leg
- Easy to use



## AML Pantyhose Aid

### SLDA5305

- Lightweight and easy to put on
- Flexible plastic trough can be bent as necessary to fit any size stockings
- Slide feet through and pull the two straps up until the pantyhose can be pulled up within reach




## Aluminium Reacher

### MP30801 (660 mm)

### MP30802 (815 mm)

- A great gadget to help pick up objects from the floor or out of reach
- Lightweight anodized aluminium body
- Plastic handle
- Magnet in jaw for picking up light metal objects



## Nipper Reacher

### MP30803 (600 mm)

### MP30804 (750 mm)

- Multi-Use Reacher with a lightweight anodized aluminium frame and plastic handle
- Magnet in jaw for picking up light metal objects and hook on upper jaw retrieves articles of clothing
- Slip resistant surface on jaws for secure hold
- Jaws rotate to avoid excessive wrist movement





## Adjustable Back Rest

### MP30101

- A supportive backrest designed to make sitting in bed easier
- Made of a white powder coated frame with a breathable nylon mesh fabric
- Comes complete with a padded, removable head rest
- Has five adjustments of angle which allow it to go from a low lying position to an upright sitting position
- For storage, it folds flat



## AML Coccyx Orthopaedic Cushion

### AMLCOC


- Coccyx orthopaedic cushion relieves tailbone pressure and encourages a correct sitting position
- Ergonomically designed from PU foam
- Unique cut section for coccyx
- Removable 100% polyester cover
- Easy maintenance



## AML Tablet Crusher

### SLHA4132

- A grinding mill for medicines.
- Turn the screw top to crush the tablets into powder
- Spare tablets can be stored in a compartment in the lid
- Suitable for carrying in a pocket and easy to clean
- Dimensions are 55 x 44 mm



## Long Plastic Shoe Horn

### PA16D019

- Strong flexible and latex free, this plastic shoe horn comes with a handy hang up hole
- Length is 18"



## Rollator - Petite 6" Wheel

### MP10103G

- A smaller rollator suitable for people 5'2" or shorter
- The cable loop brakes are easy to use
- They also have a locking feature which allows the brakes to be "locked on", necessary when using the wheeled walker's built-in seat
- A small bag under the seat allows for personal items to be safely and discreetly carried
- Maximum User Weight 136 kg



## Rollator - Bariatric 7.5" Wheel

### MP10200S

- A higher weight capacity than a regular rollator, it offers stability and mobility for the larger adult
- The rustproof, robust aluminium frame is extremely durable
- The heavy duty castors can be used both indoors and outdoors
- It has a soft padded seat and a removable backrest
- Folds up for ease of transport
- Some assembly required
- Maximum user weight 226 kg



## AML Gutter Frame Rollator

### MP61040070

- Forearm Rollator, which is great for people with weaker arms and shoulders
- Black and grey steel frame
- Folding walker designed to support the user in an upright position with padded armrests for comfort.
- Features a narrow seat and cable brakes on the rear wheels
- Gutter-style arm supports are height adjustable.
- Tray that can be slotted onto the seat for carrying small items
- Removable basket



## AML Four Wheeled Rollator - 6" Castors

### MP10103G

- A smaller rollator suitable for people 5'2" or shorter
- The cable loop brakes are easy to use
- They also have a locking feature which allows the brakes to be "locked on", necessary when using the wheeled walker's built-in seat
- A small bag under the seat allows for personal items to be safely and discreetly carried
- Maximum User Weight 136 kg



## AML Four Wheeled Rollator - 7.3" Castors

### AMLCOC

- Coccyx orthopaedic cushion relieves tailbone pressure and encourages a correct sitting position
- Ergonomically designed from PU foam
- Unique cut section for coccyx
- Removable 100% polyester cover
- Easy maintenance



## Crook Handle Cane

### MP10601K

- Made of aluminium, this cane is lightweight and easy to use
- With a padded foam grip, it is comfortable to hold and its design makes it easy to hook onto a table or chair when not in use
- Black only
- Maximum User Weight 125 kg



## Swan Handled Cane

### MP10602B (Bronze)

### MP10602K (Black)

- This offset cane is made out of lightweight aluminium, and is easily height adjustable and features a soft foam handle
- The swan neck design helps to transfer the weight of the user, making it more stable for many people
- Maximum User Weight 125 kg



## Fischer Handle Cane


### MP10605B (Bronze Left)

### MP10605K (Black Left)

### MP10606B (Bronze Right)

### MP10606K (Black Right)

- Lightweight cane designed with arthritic grip
- The shaft is made of lightweight aluminium and it is easily height adjustable



## T-Handled Cane

### MP10603K (Black)

### MP10603M (Multi)

- Lightweight aluminium, height adjustable straight cane which features an ergonomic T shaped plastic handle
- Maximum user weight 125 kg



## Floral Foldable Cane

**MP10604V (violet)**

**MP10604P (peony)**

**MP10604L (lilac)**

- Handy foldable cane ready for when you need it
- The aluminium shaft allows this cane to be lightweight, durable and adjustable
- Comes complete with its own pouch for easy storage




## Foldable Cane

**MP10604B (Bronze)**

**MP10604K (Black)**

- Handy foldable T-handle cane ready for when you need it
- The aluminium shaft allows this cane to be lightweight, durable and adjustable
- Maximum user weight 125 kg



## AML Cane with Assist Handle

**AML11599**

- Lightweight cane with additional assist handle
- Height adjustable with pin and hole shaft
- Maximum user weight 125 kg



## Quad Cane

**MP10501K (Small Base)**

**MP10502K (Large Base)**

- This cane has been adapted for increased stability; through providing a rectangular base for support
- Made of lightweight aluminium, this cane has an offset foam coated handle
- Easily height adjustable





## Folding Seat Cane

### MP30602

- Stylish folding cane seat provides you with a sturdy support cane when closed and a comfortable seat to rest on when open
- Height adjustable cane easily unfolds to give the user somewhere to rest when needed
- Made of strong aluminium, it is lightweight and durable
- Maximum user weight 90 kg



## Cane Holders

### MP10614

- This clip enables cane to suspend from table edge or push clip up cane shaft until it meets underside of table
- Makes it easier than having to reach down and pick the walking stick up from the floor
- Black only



## Cane Tip - Black

### MP60492

- Rubber cane tip
- 19 mm diameter
- Ideal for use with MP10601, MP10602, MP10603 & MP10604 canes

### Quad Cane

#### MP60912 (Small)

#### MP60913 (Large)

- Small tip is 13 mm diameter for use with MP10501K quad cane
- Large tip is 16 mm diameter for use with MP10502K quad cane



## AML Deluxe Rubber Cane Tip - Black

### AMLDRCTBLK

- Can tip for walking stick to give added stability and allows the cane to stand alone when needed
- Rubber cane tip
- Offers the user more stability than a standard single cane tip, especially in uneven terrain
- Approx 18mm diameter



## Children's Elbow Crutches

**SMBE3640 (4 - 7 Years)**

**SMBE1070 (6 - 10 Years)**

- Strong but lightweight, these aluminium crutches have double height adjustment and are finished in a cheerful buttercup yellow
- They are completely latex free - including ferrule
- Sold as a pair
- Maximum User Weight 100 kg



## Forearm Crutches

**MP10401 (Short/Medium)**

**MP10402 (Medium/Large)**

- Made of lightweight aluminium with a moulded ergonomically contoured arm cuff for comfort
- Double height adjustment with an extra-large rubber tip that provides added stability
- Sold as a pair
- Maximum User Weight 115 kg



## Forearm Crutches – Comfy Handle

**MP10403**

- Designed for long term users
- Double height adjustment crutch made from anodised aluminium
- Features a deluxe contoured "Ortho-Grip" Comfy handle that ensures the weight is spread through the palm
- Sold as a pair
- Maximum User Weight 115 kg



## Forearm Crutches – Heavy Duty

**MP10404**

- A sturdy elbow crutch for a heavier adult
- Made of lightweight aluminium with a moulded ergonomically contoured arm cuff for comfort
- Double height adjustment with an extra-large rubber tip that provides added stability
- Sold as a pair
- Maximum User Weight 160 kg





## Elbow Crutches - Adult Tall

### SMBE3099

- Sturdy crutches for taller people
- Made of lightweight aluminium with a vinyl coated ergonomically contoured arm cuff for comfort
- Double height adjustment with an extra-large rubber tip that provides added stability
- Sold as a pair
- Maximum User Weight 180 kg





## CTS Wrist Brace

<b>AML87162 (LHS XSmall)</b>	<b>AML87152 (RHS XSmall)</b>
<b>AML87163 (LHS Small)</b>	<b>AML87153 (RHS Small)</b>
<b>AML87165 (LHS Medium)</b>	<b>AML87155 (RHS Medium)</b>
<b>AML87167 (LHS Large)</b>	<b>AML87157 (RHS Large)</b>
<b>AML87168 (LHS XLarge)</b>	<b>AML87158 (RHS XLarge)</b>

- Our top selling wrist brace - it's available from X-Small to X-Large in either Left or Right side
- Cotton elastic with circumferential strap and Velcro closure
- Preformed aluminium stay
- An ideal support for Carpal Tunnel Syndrome, as well as post fracture and wrist sprains and strains



## Allard S.O.T. Resting Splint

<b>AL287101008 (LHS XXXSmall)</b>	<b>AL287102008 (RHS XXXSmall)</b>
<b>AL287101009 (LHS XXSmall)</b>	<b>AL287102009 (RHS XXSmall)</b>
<b>AL287101010 (LHS XSmall)</b>	<b>AL287102010 (RHS XSmall)</b>
<b>AL287101011 (LHS Small)</b>	<b>AL287102011 (RHS Small)</b>
<b>AL287101012 (LHS Medium)</b>	<b>AL287102012 (RHS Medium)</b>
<b>AL287101013 (LHS Large)</b>	<b>AL287102013 (RHS Large)</b>

- Smooth and lightweight
- Allows adjustment to the desired position
- Aluminium core is embedded into polyethylene foam, covered with fabric
- Offers a resting position which can deliver pain relief, relaxation to the hand and provides a good biomechanical position that may reduce the risk of flexor shortening of the wrist and fingers



## Allard Selection Thumb Rigid Support

<b>AL352021011 (LHS Small)</b>	<b>AL352022011 (RHS Small)</b>
<b>AL352021012 (LHS Medium)</b>	<b>AL352022012 (RHS Medium)</b>
<b>AL352021013 (LHS Large)</b>	<b>AL352022013 (RHS Large)</b>

- Keeps thumb in a relaxed and pain free position
- Mild to moderate arthritis in the first CMC-joint and protects the CMC-joint and the first MCP-joint
- Incorporates two elastic materials laminated together with a core of a more rigid material to offer the required stability





## Allard Selection Wrist Support For Children

<b>AL353051010 (LHS XSmall)</b>	<b>AL353052010 (RHS XSmall)</b>
<b>AL353051011 (LHS Small)</b>	<b>AL353052011 (RHS Small)</b>
<b>AL353051012 (LHS Medium)</b>	<b>AL353052012 (RHS Medium)</b>
<b>AL353051013 (LHS Large)</b>	<b>AL353052013 (RHS Large)</b>

- Low cut shape in the palm area to allow for the bending zone
- Well cut around the thumb muscles to allow full movement of the thumb
- It comes with a simple aluminium stay
- Helps guide the hand and wrist into a functional position
- Designed to reduce spasticity in hand and wrist in neurological conditions



## Allard Selection Thumb Support For Children

<b>AL353001010 (LHS XSmall)</b>	<b>AL353002010 (RHS XSmall)</b>
<b>AL353001011 (LHS Small)</b>	<b>AL353002011 (RHS Small)</b>
<b>AL353001012 (LHS Medium)</b>	<b>AL353002012 (RHS Medium)</b>
<b>AL353001013 (LHS Large)</b>	<b>AL353002013 (RHS Large)</b>

- Elastic all the way around
- Designed to keep thumb in a relaxed and pain free position
- Designed to guide the thumb and hand into a more functional position to facilitate the gripping ability



## Ossur Wrist Brace Exoform Carpal Tunnel

<b>OS517083 (LHS Small)</b>	<b>OS517073 (RHS Small)</b>
<b>OS517085 (LHS Medium)</b>	<b>OS517075 (RHS Medium)</b>
<b>OS517087- (LHS Large)</b>	<b>OS517077 (RHS Large)</b>

- Anatomically correct fit
- Pre-shaped into a neutral position for comfortable limited motion/flexibility
- Soft touch, contoured interior ensures a comfortable fit
- Excellent for more active users - especially those involved in more extreme sports
- Ideal for Carpal Tunnel Syndrome, Tendonitis and soft tissue injuries





## Functional/Night Resting Splint

**AML35360 (LHS Small)**

**AML35361 (LHS Medium)**

**AML35362 (LHS Large)**

**AML35365 (RHS Small)**

**AML35366 (RHS Medium)**

**AML35367 (RHS Large)**

- Forearm brace, ideal for use post surgical, CVA, TBI and carpal tunnel
- Aluminium stays throughout the foam, the thumb position adjustable
- Brace can be bent as required
- Thick padding over wrist ensures comfort



## Procure Universal Wrist Wrap

**PR82300**

- Elastic all the way around
- Designed to keep thumb in a relaxed and pain free position
- Designed to guide the thumb and hand into a more functional position to facilitate the gripping ability



## Scott Universal Wrist Brace

**SC1385**

- Simple design, with cotton sleeve that allows additional comfort and easy one handed application
- Designed to fit left or right wrist
- Cotton-elastic construction with heavy-duty vinyl strip protecting wrist
- 50 mm wide elastic wrist strap for greater wrist stability, along with a removable aluminium stay
- One size fits most



## Scott Unifoam Wrist Support

SC4005LR (LHS Regular)  
SC4005LXL (LHS X-Large)

SC4005RR (RHS Regular)  
SC4005RXL (RHS X-Large)

- Extremely comfortable, soft foam wrist brace
- Ideal for night resting and for Carpal Tunnel Syndrome
- Removable aluminium stay and extra padding around key joints



## Scott Elastic Wrist Brace With Palm Stay

SC4039LS (LHS Sml)  
SC4039LM (LHS Med)  
SC4039LL (LHS Lge)  
SC4039LXL (LHS X-Lge)

SC4039RS (RHS Sml)  
SC4039RM (RHS Med)  
SC4039RL (RHS Lge)  
SC4039RXL (RHS X-Lge)

- Immobilises the wrist and alleviates pain
- Ideal for carpal tunnel syndrome
- Constructed of 175 mm beige duo-tension elastic for varied tension around the palm and wrist area



## Foam Limb Restraint

PR90110

- Contact closure with double "D" ring cinch lock
- Foam Limb Restraint, Adult, Universal Size



## Procure Universal Wrist/Thumb Wrap

### PR82600

- Simply designed product that supports both the wrist and thumb in opposition
- One piece neoprene wrap what offers wrist support with thumb abduction
- Velcro fastening that allows for easy, one-handed applications and adjustment
- Reversible for left or right wrist



## LMB Finger Splints

**LMB501AA (X-Small Size 1)**    **LMB501B (Medium Size 3)**  
**LMB501A (Small Size 2)**        **LMB501C (Large Size 4)**

- Bendable wire in closed-cell polyethylene foam
- Wire allows for easy adjustment to meet individual splinting needs and to accommodate edema
- Large curved dorsal pad helps distribute pressure
- Measure finger from palmar crease to distal finger crease



## Finger Cot - 12 Pack

**PR71903 (Small)**  
**PR71905 (Medium)**  
**PR71907 (Large)**

- Lightweight
- Padded aluminium
- Easy application and removal
- Sold in 12 pack



## Staxx Finger Splint

**PR72241 (Size 1)**                    **PR72245 (Size 5)**  
**PR72242 (Size 2)**                    **PR72246 (Size 5.5)**  
**PR72243 (Size 3)**                    **PR72247 (Size 6)**  
**PR72244 (Size 4)**                    **PR72248 (Size 7)**

- Protect and support the distal joint of the finger in extension
- Permitting unrestricted movement of the proximal interphalangeal joint
- Sold individually and comes in skin tone



## Procare Contoured Wrist Support

### PR87007

- Perforated suede construction with soft flannel lining and cotton stockinette for patient comfort
- D-rings with loop/lock closure provide a secure fit
- Aluminium palmer stay can be contoured to fit each patient
- Designed to prevent brace from extending beyond palm crease
- Ideal for Carpal Tunnel Syndrome and immobilisation of the wrist



## Procare Abducted Thumb Spica

### PR87100

- Volara foam laminate construction with contact closure
- Provides support to the thumb and wrist
- Can be used for soft tissue injuries to the thumb, tendonitis, arthritis or deQuervain's Syndrome
- One size fits most



## Scott Wrist Wrap

### SC1307

- Made of breathable elastic with Velcro closures
- Wrap around strap has double locking system for desired tension and comfort
- 3" wide, beige in colour, one size fits most



## Scott Universal Tennis Elbow Support

### SC1969

- All round tennis and golfers elbow support
- 50 mm wide Royal blue plush fabric on both sides, accepts Velcro closure anywhere
- Easy to adjust with loop lock ring



## AML Elbow Ranger

**AML94210 (Medium)**

**AML94220 (Large)**

- For post-operative immobilisation and protected range of motion
- Elbow Ranger is ideal when a neutral position is desired
- For non-operative application, collateral ligament strains and ruptures and ulnar nerve transposition
- May be used following cast removal to stabilise fractures and for post-surgical immobilisation
- Removable handgrip assists in the prevention of pronation/supination
- Aluminium R.O.M. bilateral hinges allow for adjustments in 30 degree increments with support by the formed Kydex shells



## AML Shoulder Abduction Sling

**AML31201 (Small)**

**AML31202 (Medium)**

**AML31203 (Large)**

**AML31204 (X-Large)**

- Fits left or right
- Allows shoulder and arm to remain in a neutral position
- Pillow allows auxiliary air exchange to reduce the risk of secondary infections
- Breathable Fabric and Padding Sling for comfort



## AML Clinic Shoulder Immobiliser

**AMLCSIXSML (X-Small)**

**AMLCSISML (Small)**

**AMLCSIMED (Medium)**

**AMLCSILGE (Large)**

**AMLCSIXLG (X-Large)**

- Simple to use, multipurpose shoulder immobiliser
- Sized canvas sling and soft foam waist strap immobilise the shoulder while supporting the elbow, wrist and hand





## SAE - Spinal Anterior Extension

### AMLSAE42215

- Thoraco Lumbar Sacral Orthoses (TLSO)
- Comfortable, lightweight, universal sized
- Anatomically shaped Hyper-extension Orthosis
- Low profile and very easy to fit
- Constructed of anodized aluminium. Kydex plastic pads lined with soft foam for easy modification with a heat gun
- Hinged sterna and pubic assembly
- Ideal for compression fractures, kyphosis, osteoporosis and other conditions that require hyper-extension of the spine



## Procure Clavicle Splint

### PR85003 (Small)

### PR85007 (Large)

### PR85005 (Medium)

### PR85008 (X-Large)

- Adjustable prong buckles ensure correct fit and prevent slippage
- Stockinette-covered felt shoulder straps
- Ideal for clavicular fractures and postural problems
- Washable



## AML Multifit Back Support

### AMLMULTIFITB (Black)

- Basic back support with double elastic pulls on each side
- Ideal for short-term use providing a very economic solution
- Fits waist 800 - 1125 mm



## Scott Universal 6" Sacral Support

### SC903

- Support of the lower back and spine
- 6" wide elastic and Durafoam, which combine to provide support to the sacral area
- Adds pressure to the lower spine to help stabilise the spine and relieve back pain
- 3" wide side tension strap for added compression
- Universal Size - Fits hip circumference of 700 - 1200 mm



## Scott Elastic Duo-Adjustable Back Support

**SC3245S (Small)**

**SC3245L (Large)**

**SC3245M (Medium)**

**SC3245XL (X-Large)**

- Elastic back support provides excellent abdominal and lumbosacral support
- 225 mm wide elastic back panel with 75 mm wide elastic duo-tension straps
- Ideal for women and long term use



## Scott Scrotal Support

**SCSA1503S (Small)**

**SCSA1503M (Medium)**

**SCSA1503L (Large)**

- Scrotal athletic support
- The fabric is comfortable and breathable
- Measure around waist for size
- Hand wash with a mild detergent in cold water



## Procure Pavlik Harness


**PR90512 (XSmall)**

**PR90515 (Medium)**

**PR90513 (Small)**

**PR90517 (Large)**


- For treatment of hip dysplasia and dislocation
- Chest band and shoulder straps made with heavy webbing with Kodel Polyester pile backing for durability and comfort
- Follow a qualified healthcare provider's instructions regarding usage and infection control



## Procare Overdoor Cervical Traction Kit

### PR90130

- Relieve neck pain for anyone suffering from neck arthritis, herniated/bulging disc in the neck, neck strains and cervical muscle spasms
- Kit can be used at home effectively and safely
- Includes overdoor bracket, deluxe head halter, water weight bag, spreader bar and cord to administer cervical traction at home or the clinic



## Ossur Miami J Collar Set

### OSMJ-300 (Short)

### OSMJR-300 (SML W/REPL PAD)

### OSMJR-400 (MED W/REPL PAD)


- Offers superior immobilisation that simplifies size selection by having 6 unique front and back components
- Engineered to work optimally together for various precautions



## Miami J Collar Thoracic Extension

### OSMTR900 (W/REPL PADS)

- Fits all of the Miami J collars and is ideally used whenever greater stability of the cervical junctions is required and the control of the collar is needed to be extended
- Fitting is tool free and it is metal free to allow for ease of imaging
- Please note that the collar needs to be ordered separately



## S. O. M. I. Brace

### SOMI

- Sternal Occipital Mandibular Immobilizer, designed to provide flexion, neutral or hypertension positioning and immobilization to the neck and spine
- Very adjustable and has a removable chin plate which helps the patient with their hygiene and also eating
- It is made of vinyl-coated aluminium so it is lightweight and can also be contoured to the user



## Philadelphia Adjustable Collar

### OSPHP-AA

- Immobilisation and comfort
- Anterior height adjustability accommodates neck height from 2¼ to 4¼ inches
- Tear-off tabs accommodate neck circumference from 10 to 19 inches
- Restricts cervical spine flexion, extension and rotation
- Total cervical arch support maintains neutral alignment of the c-spine
- Water-resistant for bathing and aqua-therapy



## Philadelphia CoolMax Neck Liners

### OSPHPLCS (Small)

### OSPHPLCM (Medium)

### OSPHPLCL (Large)

- Sold separately



## Philadelphia Tracheotomy Rigid Collar

### OSPHP-TXX

(available a a broad range of sizes from Infant to Adult)

- Provides strong, open, anterior support
- Constructed of laminated foam plastic
- Recommended for single patient use



## Patriot Extrication Collar

### OSPHP-PP

- Either a supine or seated position
- Restricts cervical spine flexion, extension and rotation
- Adjustable bilateral Velcro® hook and loop fasteners
- Polyethylene Front and Back pieces for improved rigidity
- Large trachea opening for quick access to carotid pulse monitoring and emergency tracheotomies



## Headmaster Collar

### SDHMC03-G (Small),

### SDHMC02-G (Medium)

### SDHMC01-G (Large)

- Improves breathing and swallowing by opening the trachea and esophagus
- Cool and comfortable brace, with low heat retention
- Provides firm support for the head
- Quick release velcro makes it easy to put on and take off



## AML Lower Leg Walker

**AMLLLWS (Small)**

**AMLLLWM (Medium)**

**AMLLLWL (Large)**

- Low profile walker that provides stability for lower leg injuries
- Simple and cost-effective solution that provides the support you need to rehabilitate your injury
- It is made from a breathable fibre liner and a non-skid rocker sole
- Comes in three sizes and is universal to fit your left or right foot



## AML R.O.M. Walker

**AMLRMWS (Small)**

**AMLRMWL (Large)**

**AMLRMWM (Medium)**

- Low profile brace with Dual Action Hinges
- Adjustable strapping to give extra support.
- Accordion moulding in the popliteal area eliminates bunching and increases comfort
- Neoprene lined with Lycra to offer a breathable contact area for comfort and extensive use



## Ossur Formfit Night Splint With Tread

**OS50023 (Small)**

**OS50025 (Medium)**

**OS50027 (Large)**

- 3-dimensional padding throughout the shell
- Adjustable flexion straps that control the amount of stretch provided for effective treatment and provide a gentle stretching of the plantar fascia and Achilles tendon
- Is ideal for treatment of Plantar Fasciitis, Achilles Tendonitis and other lower extremity overuse injuries



## Ossur Walker Equaliser

**OSW0600BLK**

- Unique shock absorbing sole minimises the impact on the heel while walking, preventing further irritation of injuries
- Rocker bottom walking base on this moonboot is specially designed to be low profile and wide to promote a natural and stable gait
- It is ideal for soft tissue injuries, grade 2 and 3 sprains, or stable fractures as well as post-operative, trauma, or rehab use





## Ossur Airform Night Splint

**OS50045 (Medium)**

**OS50047 (Large)**

- Designed to maintain the ankle in a neutral position (right angle, as if standing)
- Provides a gentle stretch on the Plantar Fascia and Achilles Tendon
- Lightweight, low profile slim design makes it comfortable to wear at night which means better compliance
- Effective for treating plantar fasciitis, achilles tendonitis, achilles strains, achilles bursitis, drop foot and post-static pain



## Ossur Formfit Ankle Brace (Speed Brace)

**OSB212000001 (X-Small)**

**OSB212000004 (Large)**

**OSB212000002 (Small)**

**OSB212000005 (X-Large)**

**OSB212000003 (Medium)**

- Provides the compression and comfort of a soft ankle support with 20 - 35% more protection than traditional ankle stirrups
- Exoform technology forms a plastic exoskeleton directly onto soft breathable fabric
- Speed laces allows the user to lace the brace up in seconds, saving time and encouraging the brace's use



## Procure Neoprene Ankle Sleeve

**PR82303 (Small)**

**PR82305 (Medium)**

**PR82307 (Large)**

- An 1/8" neoprene sleeve providing warmth and compressive support. for ankle sprains and soft tissue injuries
- The style is easy-on and off



## Allard ToeOFF 2.0

**AL289221011 (LHS Small)**

**AL289221012 (LHS Medium)**

**AL289221013 (LHS Large)**

**AL289222011 (RHS Small)**

**AL289222012 (RHS Medium)**

**AL289222013 (RHS Large)**

- Carbon composite dynamic Ankle Foot Orthosis solution designed for users with mild to severe foot drop
- Organic shape, with short side extensions and more rounded shape
- Carbon and Kevlar design absorbs energy at the heel strike and returns it at toe-off, restoring a more normal and dynamic gait



## Allard Blue Rocker 2.0

**AL289401012 (LHS Medium)**

**AL289401013 (LHS Large)**

**AL289402011 (RHS Small)**

**AL289402012 (RHS Medium)**

**AL289402013 (RHS Large)**

**AL289402014 (RHS X-Large)**

- Carbon composite dynamic Ankle Foot Orthosis solution designed for users with mild to severe foot drop
- Organic shape, with short side extensions and more rounded shape
- Carbon and Kevlar design absorbs energy at the heel strike and returns it at toe-off, restoring a more normal and dynamic gait



## Allard SoftKIT 2.0

**AL287500011 (Small)**

**AL287500012 (Medium)**

**AL287500013 (Large)**

- For ToeOFF 2.0 and Blue Rocker 2.0 AFO
- Soft interface between the orthosis and the tibial crest
- Padding system that includes two narrow neoprene pads to create a pressure-free channel
- Helps relieve any potential discomfort from bony prominences



## Knee Immobiliser

**PR80023 (Small)**

**PR80025 (Medium)**

**PR80027 (Large)**

**PR80028 (X-Large)**

- Comfortable, sturdy, foam/mesh construction with alternating web straps and loop/lock contact closure
- All sizes are 20" / 50.8 cm in length
- Comfortable and sturdy
- Loop/lock contact closure for easy fitting
- Combination of stays for proper immobilisation
- HI-beam stays
- Foam/mesh construction with alternating web straps



## Ossur Wraparound Hinged Knee Support

**OS212003 (Small)**

**OS212008 (X-Large)**

**OS212005 (Medium)**

**OS212009 (XX-Large)**

**OS212007 (Large)**

**OS212010 (XXX-Large)**

- Low profile brace with Dual Action Hinges
- Adjustable strapping to give extra support
- Accordion moulding in the popliteal area eliminates bunching and increases comfort
- Neoprene lined with Lycra to offer a breathable contact area for comfort and extensive use



## Knee Ranger

**PR94240**

- Ideal for providing protected range of motion for postoperative immobilization
- The dual-axis, polycentric hinge allows for convenient dial-in adjustment with flexion and extension stops at 0°, 15°, 30°, 60° and 90°.
- Trimmable foam leg wraps for a customised fit
- Medial/lateral plastic shells riveted to hinge bars provide additional support following post-operative procedures
- Accommodates Ankle/Shoe Insert for added rotational control



## Lace Up Ankle Brace

**PR81313 (Small)**

**PR81317 (Large)**

**PR81315 (Medium)**

**PR81318 (X-Large)**

- Fits easily into any athletic shoe
- Sturdy canvas stirrup brace
- Designed with removable plastic medial / lateral stays
- Provide support for acute ankle ligament sprains and chronic instabilities



## Scott Double Strap Ankle Support

**SC325S (Small)**

**SC325L (Large)**

**SC325M (Medium)**

**SC325XL (X-Large)**

- Figure of 8 Ankle Brace
- Slip-on ankle support with open heel, made of beige four-way stretch spandex
- Double strap style offers compression to both medial and lateral sides of the ankle
- Provides mild ankle protection to help prevent injury or re-injury for chronic ankle instabilities or strains
- May help control over rotation of the ankle.



## AML Supralite Thermoplastic Ankle Foot Orthosis (A.F.O.)

**AMLSUPAFOFL (Female LHS)**

**AMLSUPAFOML (Male LHS)**

**AMLSUPAFOFR (Female RHS)**

**AMLSUPAFOMR (Male RHS)**

- Made of durable, orthopaedic grade thermoplastic
- Can be custom trimmed and reshaped using cast scissors and a heat gun
- This product is recommended for single patient use



## Black Dictus Band - Outdoor

**ERDS02J (Junior)**

**ERDS02M (Medium)**

**ERDS02L (Large)**

**ERDS02XL (X-Large)**

- An orthosis to help people who suffer from Drop Foot
- Drop Foot (or Foot Drop) is a neuromuscular (nerve and muscle) disorder that affects the patient's ability to raise their foot at the ankle
- Included parts: 1 x Ankle band, 2 x Latex Rubber Bands,
- 1 x Monohook, 2 x Steel hooks



## Replacement Mono Hook for Dictus Band

**ERDS0102**

- 'Hook' system that can be used with the outdoor Dictus Band
- It slides into the tongue area of the shoe underneath the laces, with the hook coming through the laces
- Latex Rubber Band then hooks into the hook of the Mono Hook and also the hook on the Dictus Band to assist with drop foot
- Sold individually



## Replacement Spring Steel Hooks for Dictus Band - Pair

**ERDS0103**

- For use with an outdoor Dictus Band to hook through shoe laces and attach Latex Rubber Bands to
- Sold as a pair



## Replacement Latex Rubber Band

**ERDS01041**

- For use with Dictus Bands - either indoor or outdoor
- Fantastic to have on hand for when your current bands wear out over time
- Sold individually





## Black Dictus Band Shoeless Cuff

**ERDS03L (LHS)**

**ERDS03R (RHS)**

- Simple AFO for Drop Foot that can be used without shoes in conjunction with the ankle strap from the Outdoor Dictus Band or Dictus Flex
- Please note: Only lower foot cuff included, Black Dictus Band (sold separately)



## Tulis Classic Gel Heel Cups

**TULIS10215 (Regular)**

**TULIS10216 (Large)**

- Tulis Gel is a specially developed polymer that is lighter, softer and more resilient than the material in similar heel cups
- The patented waffle design provide the ultimate shock absorption for use in shoes of all types
- Sold as a pair



## Tulis Heavy-Duty Gel Heel Cups

**TULIS10219 (Small)**

**TULIS10217 (Regular)**

**TULIS10218 (Large)**

- Tulis Heavy-Duty Gel heel cups provide the ultimate in shock absorption for high-impact activities
- Sold as a pair



## Reciprocating Pulley Kit

**ERDS03L (LHS)**

**ERDS03R (RHS)**

- Simple AFO for Drop Foot that can be used without shoes in conjunction with the ankle strap from the Outdoor Dictus Band or Dictus Flex
- Please note: Only lower foot cuff included, Black Dictus Band (sold separately)



## Pedal Exerciser

**MP30401**

- Created especially for people with limited mobility, who cannot use conventional exercise equipment
- The pedal exerciser helps to strengthen your legs and/or arms and hands while increasing your blood circulation
- Chrome plated with full tension range



## Gemstim TENs Machine

### GEMSTIM

- A cost effective TENs unit ideal for single patient use
- The unit is a dual channel electric stimulator for active treatment applications
- Comes complete with: Four pre-wired reusable self adhesive electrodes; 9-volt battery; Leads (two channel); Instruction manual; compact carry case



## Self Adhesive Electrodes TENS machine

### ELECTSAD

- Sold separately



## Obstetric TENS Machine

### GEMSTIM

- A cost effective TENs unit ideal for single patient use
- The unit is a dual channel electric stimulator for active treatment applications
- Comes complete with: Four pre-wired reusable self adhesive electrodes; 9-volt battery; Leads (two channel); Instruction manual; compact carry case



## Self Adhesive Electrodes Obstetric TENS machine

### OBSELECT

- Sold separately



## Leads For Gemstim TENs Machine

### LEADSGS

- Replacement electrode leads to suit GEMSTIM TENS machines
- Sold as a pair



## Gel for TENs Machines

### GEL-60ML

- Gel for use with TENS Machines that use silicone electrodes
- Can also be used with products such as Foetal Dopplers
- 60 ml tube



## McDavid Reusable Cold or Hot Pack

### MCD DAVID ICE

- Low cost, reusable ice pack or hot compress
- Filled with non-toxic gel that can be boiled or microwaved as well as frozen
- Clear instructions printed on pack for heating or freezing
- Ideal for clinic use



Please visit our website  
[www.alliedmedical.co.nz](http://www.alliedmedical.co.nz)  
for our full range of products



**0800 31 61 81**  
**helpis@alliedmedical.co.nz**

July 2023. This information is correct at the time of printing.  
From time to time, product specifications are changed by manufacturers without prior notice.



**ALLIED  
MEDICAL  
LIMITED**