

Patron Piper Comfort Buggy

The Piper Comfort buggy is used for short term transportation and activities for your child. It has a fully foldable undercarriage made from lightweight aluminium profiles.

The buggy is equipped with a seat base that has an adjustable calf rest, extended reclining backrest, solid wheels, detachable and adjustable footrest and a 5-point harness with crotch belt and shoulder padding. The variability of this buggy and the ability to attach a wide range of accessories allows the user's varying needs to be safely catered for.



Standard Features:

- Adjustable recline on the backrest
- Height adjustable footrest that can be detached
- 5-point immobilisation harness with crotch belt and shoulder padding
- Detachable wheels and forks
- One-piece foldable handlebar
- Load capacity up to 50kg



Seat Positioning



Footrest Positioning



Backrest Recline



**ALLIED
MEDICAL
LIMITED**

Specifications

Backrest Height	68cm
Seat Depth	29cm
Seat Width	35cm
Seat to Footrest	25-37.5cm
Backrest Angle	90° - 135°
Seat Angle	26°
Total Length	108cm
Total Width	59cm
Weight	14.5kg
Maximum Load	50kg

Accessories



Flexipads Set



Adjustable Headrest



Abduction Block



Feet Fixation



Classic Canopy



Footmuff - Summer



Footmuff - Winter



Raincover



Seat Inlay