

DISCOVER THE

CATALYST *ADVANTAGE*



CATALYST 5




CATALYST


Since 2006, Catalyst has been the best performing and most durable ultra-lightweight folding chair on the market. Catalyst continues to provide one-of-a-kind features, enhancements and exclusive innovations to the complex rehab industry.

Catalyst is designed to make your team more efficient. Every Catalyst model is constructed using the same high-quality components and innovative design.

Folded Width - 4"




BETTER *OPTIONS & ACCESSORIES*




SUPER LOW SWING AWAY

Lower the client's seat by 1.5 inches, making this the perfect option for foot propulsion.



HANGER RELEASE CHOICES

Three options to meet the client's or caregivers needs. Choose from a classic hanger, 4-way inline or 4-way out.



X-HINGE DUAL CROSS BRACE

The dual interlocking cross brace provides the rigidity needed to maximize the client's efficiency.



MAXX PERFORMANCE

The superior rigidity, performance and responsiveness gives the client more out of every push.




PRO ELR

Enhanced adjustability in the pivot point gives the client a personalized fit that can improve clinical outcomes.



RESIDUAL LIMB SUPPORT

The ease of adjustment and swing in swing out design allows the client to achieve the optimal angle and position with ease and without tools.



PERFORMANCE CASTER HOUSING

Easily adjust the client's caster housing to achieve their optimal setup with an additional 1-1.5" of seat slope.



TRANSIT WHEELS

Offer the client the ability to navigate narrow spaces and doorways with rear frame integrated Transit Wheels.

AXIOM SEATING

Experience Axiom's vast offering of backrests, cushions and custom seating that can be installed on Catalyst. Get simplified ordering through the ACS App.



CATALYST *ADVANTAGE*

K0005 Comparables*	Ki Mobility Catalyst 5	Sunrise Medical Quickie 2	Sunrise Medical QXi	TiLite Aero X	Motion Composites Helio A7	Motion Composites Helio A6
Standard Weight Capacity	136.1 kg	136.1 kg	136.1 kg	120.2 kg	120.2 kg	120.2 kg
Standard Weight †	10.4 kg	11.5 kg	13.1 kg	11.4 kg	Not Published	Not Published
Transport Weight §	5.8 kg	6.2 kg	8 kg	8.2 kg	6.2 kg	5.9 kg
Cross Brace Design	Dual w/ Power X-Hinge	Dual	3-Arm	3-Arm	Single	Single
Lifetime Warranty on Footrest	•					
Frame Material	7000 Series Aluminum	Aluminum	Aluminum	Aluminum	7000 Series Aluminum	6061 Series Aluminum
Titanium Frame Upgrade	•			•		
WC-19 Ansi-Resna Transit Approved	•	•	•	•	•	•
Make It Right Program	•	•	•			
Free Growth Program	•	•				
Heavy Duty Option (350 lbs)	•	•		•	•	•

* All information published on manufacturers website as of 2/18/2021

† 16x16 standard configuration, including rear wheels and footrests

§ 16x16 standard configuration, without rear wheels and footrests

