

PHILADELPHIA® TRACHEOTOMY COLLAR

INDICATIONS FOR USE

- Post-trauma
- Herniated cervical disc
- Cervical spondylitis
- Post-operative
- Cervical spondylosis
- Motor Neurone Disease
- RA and OA cervical spine

FEATURES AND BENEFITS

- Not made with natural rubber latex, non-toxic and hypoallergenic Plastazote® foam reduces skin irritation
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy.



REQUIRED MEASUREMENTS

- Neck circumference and height from tip of chin to sternal notch

ORDERING INFORMATION

Size	Neck circumference	6 cm / 2.25" height	8 cm / 3.25" height	11 cm / 4.25" height	13 cm / 5.25" height
X-Small	up to 10" (< 25 cm)	PHP-T2XS			
Small	10 – 13" (25 – 33 cm)	PHP-T2S	PHP-T3S	PHP-T4S	PHP-T5S
Medium	13-16" (33 – 41 cm)	PHP-T2M	PHP-T3M	PHP-T4M	PHP-T5M
Large	16-19" (41 – 48 cm)	PHP-T2L	PHP-T3L	PHP-T4L	PHP-T5L
X-Large	19" up (48 cm)	PHP-T2XL	PHP-T3XL	PHP-T4XL	PHP-T5XL

Size	Neck circumference	Height	Part#
Infant	6-8" (15 – 20 cm)	1.25" (3 cm)	PHP-TI
Pediatric	8-11" (20 – 28 cm)	1.75" (4.5 cm)	PHP-TP

ACCESSORIES

PHILADELPHIA COOLMAX COTTON LINER

Size	Neck circumference	Part#
Infant	6-8" (15 – 20 cm)	PHP-LCI
X-Small	up to 10" (< 25 cm)	PHP-LCXS
Small	10-13" (25 – 33 cm)	PHP-LCS
Medium	13-16" (33 – 41 cm)	PHP-LCM
Large	16-19" (41 – 48 cm)	PHP-LCL
X-Large	19" up (> 48 cm)	PHP-LCX

PHILADELPHIA TRACHEOTOMY ATTACHMENTS

Size	Stabilizer	Halo
Pediatric	PHP-OSP	PHP-OHP
Adult	PHP-OS	PHP-OH