

Vission Shoulder orthosis IMMO

Vission shoulder immobilizer is intended to immobilize and contain the upper limb, elbow, forearm flexed at 90 or 130 degrees as required. Especially suited for shoulder and pectoral girdle lesions.

The product is reversible and can be used either on the left or on the right shoulder as required.

Intended use:

Shoulder, arm, elbow and pectoral girdle trauma (fractures, dislocations and sprains). Shoulder and pectoral girdle surgery. Neurology: paralysis of brachial-plexus.

Item no.	Model	Size	Chest circumf.	Length upper arm
46022 0010	Shoulder IMMO	X-Small	68 - 78 cm	< 30 cm
46022 0011	Shoulder IMMO	Small	80 - 90 cm	< 38 cm
46022 0012	Shoulder IMMO	Medium	90 - 100 cm	< 38 cm
46022 0013	Shoulder IMMO	Large	100 - 113 cm	< 38 cm
46022 0014	Shoulder IMMO	X-Large	113 - 120 cm	< 38 cm
46022 0015	Shoulder IMMO	XX-Large	120 - 140 cm	< 38 cm

