

Karma

FLEXX JUNIOR

Specifically designed for functionality with multiple adjustments to suit children's needs as they grow.



ALLIED
MEDICAL
LIMITED



+ FLEXIBLE, ROBUST AND STYLISH

The Karma Flexx Junior is designed with plenty of adjustments to create flexible sitting to fit individual needs. With ellipse-shaped frame design, the Flexx is not only stylish but also robust and durable.

The configurable options of the Flexx Junior include variable seat depths, seat to floor heights, backrest heights and backrest angles, castor fork angles and centre of gravity. It can easily be set up to provide optimal comfort and performance.

The bright colour of key adjustment points are self-explanatory; users can easily operate the following functions independently: swing-away footrests, flip-back and adjustable armrests.

+ FEATURES

- Configurable dimensions
- Height adjustable flip-back armrests
- Swing-away 90 degree footrests with flip-up footplates
- Back angle adjustment (-4° to +16°)
- Padded height and angle adjustable backrest
- Reinforced double crossbrace
- Quick release pneumatic rear wheels
- Push to Lock Brakes with brake extension
- 2 point lap belt
- Heel loops
- Available as a transit wheelchair



+ ACCESSORIES

- Tray
- Height adjustable Push Handle



+ SPECIFICATIONS



Seat widths
10" (25.5 cm) to
14" (35.5 cm)




Seat depth
11" (28 cm) to
15" (38 cm)



Wheels
5" (front)
20" (rear)



Seat to floor height
14" - 22"
(35.5 - 55.9cm)



Overall width
18" - 22"
(45.4 - 55.6cm)



Folded width
12" (30cm)



Back height
14.5" (37cm)



Overall length
35" (89.6cm)



Overall height
32" (81cm)



Total weight
15.7 kg



User weight limit
75 kg

