

A Better Wheelchair, A Better Fit

# Karma KP25.2

The KP25.2 is suitable for medium and short distance travel. Lightweight and compact, it features a foldable frame and a removable compact battery pack. Perfect for your daily use, the KP25.2 is both convenient and portable.





# Why choose the KP25.2?

- Lightweight aircraft alloy frame with reinforced cross brace
- Easy to fold and disassemble (will fit into most cars)
- Interchangeable controller easily swapped to the left or right side
- Length adjustable 70 degree swing away legrests
- Swing away, height adjustable armrests
- Durable and easy-to-clean upholstery
- On a full battery charge, the KP25.2 can drive up to a 30km distance
- Programmable electronics

## Standard Features

- Flip-up, height adjustable armrests
- Transportation tie-down points (clearly marked on frame)
- VR2 60Ah Controller
- Positioning belt
- 70 degree front legrest hangers (with calf strap and wide profile castor)
- Reinforced cross braces
- Manual brake (pull to lock)
- Batteries in individual boxes
- Freewheel lever (allows for manual movement of wheelchair)

### **Specifications**



Seat widths 16", 18", 20"



Seat depth 16", 18"



Rear wheels 12"pneumatic



Front castor 8"x 2"solid



Seat to floor height 19" (optional 21")



Anti-tippers Removable



Back height
15.5"
(options: 18.5" with upholstery
extension 22" with extended
back canes)



Armrest Height 8.5" min 11.5" max



Overall width Seat width +7"



Overall Length 29" (40"with 70° hangers)



Overall height 36.5" (options: 39" with 22" back canes)



Batteries 12V 38AMP sealed gel



User weight limit 135ka



Speed Max 6km/h



## Compact and Manoeuvrable

One of the daily difficulties powerchair users encounter is manoeuvring in and through tight spaces like doorways and elevators. Most standard powerchairs are 25" wide and given that some doorways are as narrow as 28" or smaller, it's a tight squeeze getting through those doorways. The KP-25.2 is a compact medium-size powerchair with an overall width of up to 27", meaning it is easier to manoeuvre around tight corners and through narrow doorways.