



Karma KP25 Powerchair

The KP25 powerchair offers numerous options and accessories creating independence for the user. Its simplicity and design makes it a breeze to use!



***Lifetime
Warranty
on Frame***

***Standard
Rear Wheel
Drive
Powerchair***



**ALLIED
MEDICAL
LIMITED**

Designed especially for New Zealanders



+ 8 REASONS WHY THE KP25 IS A GREAT OPTION

1. Lightweight aircraft alloy frame with reinforced cross brace
2. Easy to fold and disassemble (will fit in to most cars)
3. Interchangeable controller – easily swapped to the left or right side
4. Length adjustable 70 degree swing away legrests
5. Swing away, height adjustable armrests
6. Durable and easy-to-clean upholstery
7. On a full battery charge can drive up to a 30km distance
8. Programmable electronics

+ SPECIFICATIONS



Seat width

16", 18", 20"



Seat depth

16", 18"
Options: 1" and 2" (+/-)
customisation for
16"-20" width frames



Rear wheel

12" pneumatic



Front castor

8" x 2" solid



Seat to floor height

19"
Options: 21"



Anti-tippers

Removable



0800 31 61 81

+ WHAT CAN'T YOU DO WITH A KARMA?

It's hard to believe, but we can assure you it's true! There is not much you can't do with this highly adjustable KP25 Powerchair. It has a comprehensive range of accessories available that suit a broad range of configurations and user needs.

Standard features



Flip back, height adjustable desk armrests



Transportation tie-down points (clearly marked on frame)



Shark controller (can be easily swapped from left to right)



Positioning belt



70 degree front legrest hangers (with calf strap and wide profile castor)



Reinforced cross braces



Manual brake (pull to lock)



Batteries in individual boxes



Freewheel lever (allows for manual movement of wheelchair)



Foldable frame (easily transportable)

Don't see what you need? Just call the Allied Medical team as we're bound to have more great ideas.



Back height

15.5"
Options: 18.5" with upholstery extension
22" with extended back canes



Armrest height

8.5" min
11.5" max



Overall width

Seat width + 7"



Overall length

29" (40" with 70° hangers)



Overall height

36.5"
Options: 39" with 22" back canes



Batteries

12V 38AMP sealed gel



User weight limit

115 kg



Speed

Maximum of 6kms per hour

+ OFFERING YOU THE COMPLETE PACKAGE

We realise that in order to find the total solution for your clients, sourcing the right powerchair seating options and accessories is also important. So we've put together a complete range of these products from some of the world's leading suppliers so you can get it all from one supplier.

ACCESSORIES AVAILABLE FROM ALLIED MEDICAL	 BACKS	 CUSHIONS	 POSITIONING BELTS	 TRAYS	 DROPSEAT BASES
 BODY POINT			✓		
 QUALITY FOR LIFE	✓	✓			✓
 parsons adl inc. <small>accessories for daily living</small>			✓	✓	
 VARILITE	✓	✓	✓		✓

For more great solutions and all the right answers, call Allied Medical's Powerchair Category Specialist on **0800 31 61 81**.

It's Wheelie Easy
with *Allied Medical*



www.alliedmedical.co.nz



Allied Medical Ltd
PO Box 302 250,
North Harbour, Auckland

Phone: 0800 31 61 81

Fax: 09 415 1686

Email: greatwheelchairs@alliedmedical.co.nz