Kimba Family





Products for combinations



Kimba seating unit



Kimba outdoor mobility base



Kimba Cross



Kimba Inline



Kimba Home



Kimba indoor mobility base

What	What do you want to order?						
	Kimba seating unit						
	Kimba outdoor mobility base						
	Kimba Cross						
	Kimba Inline						
	Kimba Home						
	Kimba indoor mobility base						
	Robby bathing support						

	Dimensions kid	Dimensions Kimba size 1 Dimensions Kimba					
① Back length [mm]		410 - 610	560 - 700				
② Seat depth [mm]		190 - 310 (140 shortened)	240 - 400				
Seat width [mm]		200 - 300	260 - 350				
③ Lower leg length [mm]		190 - 310 (122 shortened)	200 - 370				
Load capacity		Seating u Outdoor mobil					
a) User weight							
b) Medical Equipment							
* With option "footrest angle adjustable for short lower leg length" the maximum load decreases to 35 kg.							

Kimba Seating unit Order no.:

□ Provide estimate only □ Order □ Attachment for custom fabrication enclosed Billing address/customer number Delivery address Sales person Name Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/ service manuals.

Kimba Seating unit

Maximum load capacity: 40 kg
Article number: 470G71=5_AA02_C



Figure includes accessories

Standard Equipment

- · Seat and back support
- Seat width, seat depth, back support height and lower leg length are continuously adjustable
- Seat inclination and knee angle continuously adjustable
- Back support angle adjustable
- Optimized quick-release mechanism to remove the seat
- Position the seat in/against the direction of travel
- · Leg support assembly with flip-up foot plate
- Lateral side supports
- Adjustable pelvic laterals
- · Frame colour: black
- · Padding colour grey

Notice:

This product can be combined with outdoor mobility bases Kimba 470G71=5_AA03_C, Kimba Cross and Kimba Inline and also with indoor mobility bases Kimba Home and Kimba. Please fill out the corresponding order sheets separately.

Size		
KC01	Size 1 (from approx. 1 to 6 years)	-
KC02	Size 2 (from approx. 4 to 10 years)	Add. charge

-
dd. charge
_
dd. charge
dd. charge
dd. charge
dd. charge
dd. charge

Pelvic belts					
KP02B/KP02		Lap belt, padded	□ enclosed	mounted	Add. charge
KP18B/KP18		4-Point paediactric lap belt, with plastic buckles (width: 2,5 cm; length: 100 cm)	□ enclosed	□ mounted	Add. charge
KP06B/KP06		Groin strap	□ enclosed	 mounted	Add. charge
KP15B/KP15		Paediatric lap belt, with buckles (width: 2,5 cm; length: 100 cm)	enclosed	mounted	Add. charge
Chest and shoulde	er supports				
KP16B/KP16	VE	Shoulder harness Size S	enclosed	☐ mounted	Add. charge
KP17B/KP17		Shoulder harness Size M	☐ enclosed	□ mounted	Add. charge
KP19B/KP19	(B)	Chest and shoulder pad Size S	☐ enclosed	☐ mounted	Add. charge
KP20B/KP20	124	Chest and shoulder pad Size M	☐ enclosed	☐ mounted	Add. charge
KP21B/KP21		Slimline sternum harness	□ enclosed	□ mounted	Add. charge
KP22B/KP22		Dynamic harness Size XS	□ enclosed	☐ mounted	Add. charge
KP23B/KP23		Dynamic harness Size S	☐ enclosed	□ mounted	Add. charge
KP24B/KP24	1	Dynamic harness Size M	☐ enclosed	□ mounted	Add. charge
KP03B/KP03		4-point belt, padded	□ enclosed	 mounted	Add. charge
KP05B/KP05		5-Point belt, small buckle (recommended for seating unit size 1)	enclosed	□ mounted	Add. charge
KP04B/KP04		5-Point belt, large buckle	☐ enclosed	□ mounted	Add. charge
KP10B/KP10		Positioning vest	□ enclosed	□ mounted	Add. charge

Preassembly									
KP14			Mounting kit strap/vest, mounted	Add. charge					
Trunk support									
KD04B/KD04			Thoracic support	Add. charge					
KD31B			Swing away thoracic supports (laterals), enclosed	Add. charge					
KD10B			Lumbar support padding, enc- losed	Add. charge					
KD11B		4	Neck support, enclosed	Add. charge					
KC10		8	Seat depth reducer (seat depth possible up to 14 cm)	Add. charge					
KC13B			Inlay small, enclosed (spacer fabric grey, only for size 1)	Add. charge					
Arm supports	;								
KE01		F	Arm support, Standard, mounted	Add. charge					
KE02			Arm support hight/low, adjustable,mounted (not possible with country of destination Germany)	Add. charge					
Tray									
KS01B			Tray, grey, enclosed (only possible in direction of travel and in combination with KE 02 "arm support high/ low" in lowest setting only)	Add. charge					
KS02B			Tray, transparent, enclosed (only possible in direction of travel and in combination with KE 02 "arm support high/ low" in lowest setting only))	Add. charge					

Leg and foot	Leg and foot support							
KB03			Leg support, Standard, mounted (from LLL 19 cm)			-		
KB05		Page 1	Leg support, angle-adjustable, mounted (from LLL 19 cm)			Add. charge		
KB10			Leg support, angle-adjustable for short lower leg length, mounted (only for size 1, from LLL 14 cm)			Add. charge		
Accessories I	eg ar	nd foot support						
KB08B/KB08			Slip resistant footrest cover (not mounted in combination with KB10)	□ enclosed	□ mounted	Add. charge		
КВ09В			Padding, enclosed			Add. charge		
KB12B		90	Foot support – sandals, enclosed			Add. charge		
KB07B		00	Velcro footstraps, instep padded			Add. charge		
KC06B/KC06		TED	Abductor size 1, only for seat size 1, (L 9 cm × W 6 cm × H 9 cm)	 enclosed	☐ mounted	Add. charge		
KC07B/KC07		JA ?	Abductor size 2 (L 10 cm × W 7,5 cm × H 10 cm)	☐ enclosed	☐ mounted	Add. charge		
KC08B/KC08			Abductor size 3, only for seat size 2 (L 12 cm × W 9 cm × H 10 cm)	☐ enclosed	□ mounted	Add. charge		
Canopy								
KW13			Canopy in grey-melange, inside grey, incl. rain cover and mounted adapter			Add. charge		

Slip sack			
KW01B		Summer/winter slipsack Size S, enclosed	Add. charge
KW02B		Summer/winter slipsack Size M, enclosed	Add. charge
KW03B		Summer/winter slipsack Size L, enclosed	Add. charge
Accessories	seat		
KQ04B		Grip bar with padding, enclosed	Add. charge
KD20		Back angle limitation, mounted	Add. charge
KW04B		Rain cape, grey, enclosed	Add. charge
KD12B		Incontinence cover, enclosed	Add. charge
KC16B		Retrofit kit for Kimba Neo street frame (without tilt; not possible in conbination with custom fabrication platform for respirators up to 90 kg/ 198,4 lb)	Add. charge

Kimba outdoor mobility base Auftrags-Nr.:

□ Provide estimate only □ Order □ Attachment for custom fabrication enclosed Billing address/customer number Delivery address Sales person Name Date Signature

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/ service manuals.

ottobock. Kimba outdoor mobility base

Maximum load capacity: 55 kg (including seating unit)

Article number: 470G71=5_AA03_C



Standard equipment

- Folding outdoor mobility base made of high-Folding outdoor mobility base made of high-grade aluminium
- · Optimized rear axle suspension
- · Push bar with ratchet joint, height-adjustable
- Seat tilt with tilt activation at the frame incl. locking pin
- Integrated storage space (also usable as a platform for respirators)
- Automatic folding lock
- Optimized seat adapter with automatic tilt limitation
- · EVA tyres with frictionbrake
- Front wheels 7", swivelling with integrated swivel lock
- Integrated anchor points for mounting to wheelchair accessible vehicle

Size			
KA01	Size 1 (from approx.	. 1 to 6 years)	
KA05	Size 2 (from approx.	. 4 to 10 years)	
Adapter Sitzschale			
KC03		Seating shell interface	Add. charge
KB01		Leg support assistance for seating shell interface	Add. charge
KC14	4	Universal adapter (seat size 1 to base size 2) (only for country of destination Brazil)	Add. charge
KC15B		Retrofit-kit: Kimba Neo seat	Add. charge

F									
Frame colour									
013		Jet black RAL 9005 matte							
023		White RAL 9016							
Reflector sticker									
KA13			Classic	Add. charge					
KA14			Sporty	Add. charge					
Classic Sporty									
0-1									
Colour reflector sticke	er								
			Grey						
			Green						
			Turquoise						
Tyres									
KF03			EVA	-					
KF04			Pneumatic	Add. charge					
Rear wheel type									
KG01		Wheel lock (actuate	ed by the foot via wheel lock bar)	-					
KG02			anually actuated via brake levers at the push bar)	Add. charge					

Accessories						
КА07В			Frame padding, pair, enclosed	Add. charge		
KS03B			Buggy board, enclosed	Add. charge		
KA08			Tip-assist, left	Add. charge		
KA12			Tip-assist, right	Add. charge		
KS14B	Z	2	Tether band, enclosed	-		
Accessories general						
KS07B			Cupholder, enclosed	Add. charge		
KS13B			Diaper bag, grey-melange, enclosed (Contains extensive accessories such as: changing mat, dummy bag, isolation pocket for bottles and a plastic utensil bag, as well as numerous compartments and zipped pockets)	Add. charge		

| Corder no.: | Maximu base: 5tanticle | Maxim

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/ service manuals.

Name

Date

Signature

Sales person

Kimba Cross

Maximum load capacity of the outdoor mobility base: 50 kg/110 lbs (including seating unit)

Article number: HR3401=0000_K



Standard equipment

- Three-wheeled outdoor mobility base with telescoping frame, for mounting a Kimba neo seating unit, size 1 or size 2, or a custom seating shell
- Can be folded/disassembled
- Elastomer spring system
- Drum brake
- 18" rear wheels with pneumatic tyres
- Rigidly mounted 16" front wheel with pneumatic tyre, includes splash guard
- · Height-adjustable push bar
- · Frame colour: black

For children from approx. 1 to 10 years

Outdoor mobility base						
		Kimba Cross outdoor mobility base	-			
For use in conjunction with						
		Kimba seating unit, size 1				
		Kimba seating unit, size 2				
Seat adapter for Kimb	a sea	ating unit				
HR34011300		Size 1	-			
HR34021300		Size 2	-			
HR32x48300		Seating shell interface Kimba	Add. charge			

Accessories for ou	utdoor m	obility base	
HR34012100		12" swivelling wheel with quick-lock attachment	Add. charge
HR34017000		Anti-tipper, complete	Add. charge
HR34018000		Clothing guard (plug-on wheel covers)	Add. charge
HR34018100		Clothing guard extension around the wheel (not available in combination with spoke protector)	Add. charge
HR34019000		Spoke protector, transparent	Add. charge
HR34015000		Storage bag (e.g. for the swivelling front wheel)	Add. charge
HR34016100		Bicycle draw-bar, complete (includes coupling from Weber company, type E, safety strap)*	Add. charge
HR34016200		Equipment package required according to German Road Traffic Licensing Regulations (StVZO) (lighting/pennant/roll bar)	Add. charge
HR16001102		Buggy board	Add. charge
HR32xx2700-		Foot support assembly for seating shell interface Kimba	Add. charge
* According to a leafle	t for bicycl	e trailers issued by the German Federal Ministry of Transport: (www.pdeleuw.de/fahrra	d/fakten) – Recommen-

HR34016200 – Equipment package required according to German Road Traffic Licensing Regulations (StVZO)

According to a leaflet for bicycle trailers issued by the German Ministry of Transport: (www.pdeleuw.de/fahrrad/fakten)

Recommended overall weight of the bicycle trailer, including load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors. Seating shell interface and footrest assembly in combination with Squiggles and Mygo not possible.

HR34018000 – Clothing guard is a plug-on wheel cover

ded overall weight of the bicycle trailer incl. load: 40 kg. Recommended in combination with four-point belt, groin strap, spoke protectors.

HR34018100 – Clothing guard extension around the wheel cannot be used in combination with spoke protector Accessories only for corresponding size

Footrest assembly only available in combination with seating shell interface

Kimba Inline Kimba Inline Order no.: Maximum load capacity of the outdoor mobility base: 50 kg/110 lbs (including seating units) Article number: HR3255=0000_K □ Order □ Provide estimate only ☐ Attachment for custom fabrication enclosed Billing address/customer number Standard equipment Delivery address Kimba Inline twin outdoor mobility base for two Kimba seating units, size 1 Pneumatic tyres provide suspension comfort · Rear wheels with drum brake · Frame colour: black Sales person Name Twin buggy for children from approx. 1-6 years Date Signature By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/ service manuals. Inline outdoor mobility base Inline twin outdoor mobility base for up to two seating systems

positioned one behind the other, with pneumatic front wheel 200 and pneumatic rear wheel 300, complete with drum brake

Footrest assembly for front seating shell interface

Footrest assembly for rear seating shell interface

Add. charge

Add. charge

Add. charge

Add. charge

Add. charge

HR32551000-

Accessories HR16001102

HR32548300

HR32548301

HR32582700

HR32582701

Buggy board

Front seating shell interface

Rear seating shell interface

Kimba Home Order no.: □ Provide estimate only □ Order Billing address/customer number

Name

Delivery address

Sales person

Basic model

470G71=90001

471A71=SK922

ottobock. Kimba Home

Max. load capacity (including seating unit):

Size 1: 50 kg Size 2: 75 kg



- and size 2
- nit via gas djustment range
- ressure at unit adapter
- 2 with locking
- for seat angle

Add. charge

	Date Signature	 Tilt adjustment of the seat ur pressure spring 500 N (tilt a -5° till +35°) Height adjustment via gas propersion 1200 N (height of seamin. 19 cm, max. 56 cm) 4 castors (75 mm), of these device Frame colour: Kimba white Push bar incl. release lever that adjustment
nuals for accessories and products of with an Ottobock product, in addition to f the respective Ottobock products. A and accepted, our general terms and c placing an order, to conduct your own cordance with national requirements if	ply with the instructions for use/service ma- other manufacturers when combining them o the instructions for use/service manuals also you confirm to have read, understood onditions. In addition, you agree, upon conformity assessment or approval in ac- and when you implement modifications or fications of Ottobock instructions for use/	

Kimba Home indoor chassis size 1

recommended from 40 kg load

470G71=90002	Kimba Home indoor chassis size 2	
Accessories		
471A71=SK910	Anti-tip mechanism	Add. charge
471C71=SK730	Kimba seating shell interface, size 1	Add. charge
471C71=SK830	Kimba seating shell interface, size 2	Add. charge
HR32692700-000	Foot support for seating shell interface, complete, size 2	Add. charge
471A71=SK921	Gas spring (2 x 300 N) as support for seat lift function, recommended recommended from 30 kg load	Add. charge

Gas spring (2 x 500 N) as support for seat lift function, recommended

Information



Castor with locking device



Foot pedal for height adjustment



Push bar with release lever for seat angle adjustment



Seat angle adjustment



Seat heght adjustment



Adapter for Kimba seat unit



Anit-tip mechanism

Kimba indoor mobility base Order no.:

ottobock. Kimba Indoor mobility base

Article number: HR2233=0000_K





- Height-adjustable steel construction
- · Caster with swivel lock, diameter 75 mm
- · Angle-adjustable seat interface
- · Adjustable wheel base
- · Frame colour: anthracite metallic

By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.

Standard model					
HR22330000-011		Size 1	-		
HR22430000-011		Size 2	-		
For use in conjunction with					
		Kimba Neo seating unit, size 1			
		Kimba Neo seating unit, size 2			
Caster with swivel lock					
HR22321150		Diameter 75 mm	-		
HR22321160		Diameter 125 mm	Add. charge		
Seat interface					
HR32046x00		Parallel	Add. charge		
HR32x48300		Kimba	Add. charge		
Accessories					
HR22327100		Push handle	Add. charge		
HR32044x00-	\neg	Footrest assembly for seating shell interface	Add. charge		