

Patron Corzino XCountry Buggy

The Patron Corzino XCountry Buggy, a semi-adjustable rehabilitation buggy available in four sizes with a wide variety of modifications available and a high load capacity of up to 75kg. The buggy is equipped with an optimised buggy-folding system, a new padded ergonomic seat cover and a light frame with 5-point fixation ensuring safety for your child and comfortable handling. The wide variety of sizes and modifications available ensure the buggy can work well for a broad range of users.



Standard Equipment included in the base price:

- Well padded ergonomic seat cover with reinforced back
- Comfort side protection
- 2-piece 3D-Flexi footrest
- Backrest angle adjustable (2 positions)
- 5-point immobilisation harness with crotch belt and shoulder padding
- Backrest height extension (size 38, 42 only)
- 175/265 mm solid wheels
- Standard colour is red. Dark blue, green and black are also available for indent order only.

Optional equipment and accessories available at an additional cost.



**ALLIED
MEDICAL
LIMITED**

Specifications

Size	30	34	38	42
Backrest Height	54cm	67cm	85cm	88cm
Seat Depth	28cm	32cm	36cm	38cm
Seat Width	30cm	34cm	38cm	42cm
Seat to Footrest	29-42cm	30-44cm	27-43cm	33-50cm
Backrest Angle	89-90°	96-98°	94-96°	95-98°
Seat Angle	22°	23°	24°	21°
Total Length	98cm	102cm	98cm	107cm
Total Width	52cm	56cm	62cm	64cm
Weight	11.8kg	12.4kg	13.8kg	14kg
Maximum Load	45kg	55kg	75kg	75kg

Accessories



- A Headwings
- B Adjustable Headrest
- C Buggy Canopy
- D Immobilisation Waistcoat
- E H-Style Harness
- F Chest Immobilisation Belt
- G Abduction Block
- H Rain Cover
- I Two Point Safety Belt
- J Soft Work Board
- K Foldable Armrest
- L Buggy Basket
- M Footmuff (Summer or Winter Weight)
- N Extended Side Support
- O Feet Fixation



**ALLIED
MEDICAL
LIMITED**