Recliner Support System



Adapt recliners to provide secure comfort and support.

- Adjustable padded platform for proper body alignment
- Designed to safely adjust body position
- Will work with most recliner brands
- Easy to clean recliner drape and padded attachment straps protect furniture
- Includes 5-point Chest Harness, 3-point Pelvic Harness, Headrest & Leg Support



RSS folds and stores easily



The RSS unit can be used on its own or the Soft-Touch Back Liners and Hip-Flex Seat Cushions sizes 2 - 5 can be integrated into the system. These enhance lateral support and comfort. Both are sold seperately. The product can be used for both paediatric and aged care within size specifications. See over page for sizing guidelines and specs. Please note, recliner is not included.

Reclining Support System


Product Number: 66010100

RSS includes 5-point Chest Harness, 3-point Pelvic Harness, Headrest & Leg Support

Special Tomato® Recliner Support System Size Chart		
Target User Height (approx.)	35" - 74" (89-188 cm)	
Seat Width	20" (51 cm)	
Seat Depth	18" (46 cm)	
Seat Back Height	33" (84 cm)	
Seat Back Height with Seat Liners	30" (76 cm)	
Leg pad adjustment range	8" (20 cm)	
Shoulder Harness Height Levels	14" - 27" (35 - 69 cm)	
Seat Recline Range	Full range of furniture recliner	
Overall Dimensions (width x length x height)	67" x 20" (170 x 51 cm)	
Folded Dimensions (width x length x height)	36" x 21" x 10" (91 x 53 x 25 cm)	
Maximum Weight Capacity	furniture dependent	
Weight of RSS	10 lbs (4.5 kg)	

Special Tomato® Hip-Flex® Seat Cushion		
72060201	Size 2 - 10 in. d x 10 in. d - CHERRY	
72060202	Size 2 - 10 in. d x 10 in. d - LILAC	
72060204	Size 2 - 10 in. d x 10 in. d - AQUA	
72060209	Size 2 - 10 in. d x 10 in. d - DARK GRAY	
72060301	Size 3 - 12 in. w x 12 in. d - CHERRY	
72060302	Size 3 - 12 in. w x 12 in. d - LILAC	
72060304	Size 3 - 12 in. w x 12 in. d - AQUA	
72060309	Size 3 - 12 in. w x 12 in. d - DARK GRAY	
72060401	Size 4 - 14 in. w x 14 in. d - CHERRY	
72060402	Size 4 - 14 in. w × 14 in. d - LILAC	
72060404	Size 4 - 14 in. w x 14 in. d - AQUA	
72060409	Size 4 - 14 in. w x 14 in. d - DARK GRAY	
72060501	Size 5 - 16 in. w x 16 in. d - CHERRY	
72060502	Size 5 - 16 in. w x 16 in. d - LILAC	
72060504	Size 5 - 16 in. w x 16 in. d - AQUA	
72060509	Size 5 - 16 in. w x 16 in. d - DARK GRAY	

Special Tomato [®] Liner Back Cushion	
72050201	Size 2 - 8.5 in. w x 13.5 in. h - CHERRY
72050202	Size 2 - 8.5 in. w x 13.5 in. h - LILAC
72050204	Size 2 - 8.5 in. w x 13.5 in. h - AQUA
72050209	Size 2 - 8.5 in. w x 13.5 in. h - DARK GRAY
72050301	Size 3 - 10 in. w x 14.5 in. h - CHERRY
72050302	Size 3 - 10 in. w x 14.5 in. h - LILAC
72050304	Size 3 - 10 in. w x 14.5 in. h - AQUA
72050309	Size 3 - 10 in. w x 14.5 in. h - DARK GRAY
72050401	Size 4 - 13 in. w x 15.5 in. h - CHERRY
72050402	Size 4 - 13 in. w x 15.5 in. h - LILAC
72050404	Size 4 - 13 in. w x 15.5 in. h - AQUA
72050409	Size 4 - 13 in. w x 15.5 in. h - DARK GRAY
72050501	Size 5 - 14 in. w x 16 in. h - CHERRY
72050502	Size 5 - 14 in. w x 16 in. h - LILAC
72050504	Size 5 - 14 in. w x 16 in. h - AQUA
72050509	Size 5 - 14 in. w x 16 in. h - DARK GRAY



www.alliedmedical.co.nz | 0800 31 61 81 | helpis@alliedmedical.co.nz