Powerchair Specifications for Flights

All airlines have their own specific regulations, so we recommend that you contact them directly for their detailed guidelines. However, the following refers to general information that you may need to consider.

As a general rule, airlines will require all types of wheelchairs to be checked luggage. All power chairs must be stowed as check-in luggage, so it's important to minimise the possibility of damage during transit. We suggest you remove the seat cushions and any other parts that could easily become separated from the chair.

Also, it is probably a good idea to travel with some basic maintenance tools for your wheelchair.

	4FRONT	R44	Q610	J623	Q6 EDGE 2.0	Q6 EDGE Z	Q6 EDGE HD		Q6 EDGE 3
	4FRONT	R44	Q610	J023	Q6 EDGE 2.0	Q6 EDGE Z	Q6 EDGE HD	Q1450	Q6 EDGE 3
Seat widths	14", 16", 18", 20"	10" - 24" (1/2" increments)	10" – 22" (1/2" increments)	10" – 24" (1/2" increments)	10" – 24" (1/2" increments)	10" – 24" (1/2" increments)	10" – 32" (1/2" increments)	16" – 32" (1/2" increments)	12" – 24" (1/2" increments)
Seat depths	14" - 20"	10" - 24"	10" - 24"	10" - 24"	10" - 24"	10" - 24"	10" - 28"	16" - 28"	12" - 24"
Min. seat-to-floor height	19″	17″	15″	17″	15.9″	17.5″	17.5″	16″	16.5″
Overall width	Seat width + 6.5"	25″	23″	24″	24″	26"	26"	29.25"	24.4″
Overall length (Excl. footplates)	34"	43.25″	35.25″	35″	35″	36.75"	36.75"	38.25"	32.23″
Weight including batteries	97kg	105.5kg	67.4kg	97kg	97kg	113kg	112kg	135kg	80kg
Battery type	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)	Gel cell (dry)
Battery size	2 x 40 Ah	2 x 70 Ah	2 x 30 Ah	2 x 50 Ah	2 x 50Ah	2 x 70 Ah	2 x 70 Ah	2 x 70 Ah	NF-22 (55Ah)



	SKIPPI (PAED)	KP25	
Seat widths	10" - 15"	14" - 20"	
Seat depths	11.8" - 14.5"	16" - 18"	
Min. seat-to-floor height	17"	19", 21"	
Overall width	22.5"	Seat width + 6.5"	
Overall length (Excl. footplates)	27.5″	29″	
Weight including batteries	46kg	49kg	
Battery type	Gel cell (dry)	Gel cell (dry)	
Battery size	31 aH C20	2 x 40Ah	

It is important to supply your chosen airline with as much information as possible when you make your flight booking. You may need to supply them with some or all of the following information:

The type of mobility aid you are travelling with (electric or non-electric). • What type of batteries are used: gel cell (dry) or wet (spillable). • Whether the mobility aid is collapsible. • The dimensions (height, width and length) in adjusted or disassembled

- state and its total weight.

For best practise, always contact your chosen airline for their specific regulations before getting to the airport.