

## HS-520 Deluxe Mid-Range Scooter

- Maximum Speed: 12 kmph / 7.5 mph
- Capacity: 136 kg / 300 lbs
- 290mm/11" Front and Rear Wheels
- 625 mm/ 24.6" Width
- Front and rear suspension
- LED headlight
- Height and width adjustable armrest

## TECHNICAL SPECIFICATIONS

Overall Length	1200 mm / 47.2"
Overall Width	640 mm / 25.2"
Overall Height	1155 mm / 45.5"
Wheels: Front	290 mm / 11"
Wheels: Rear	290 mm / 11"
Weight w/ Batteries	94 kg / 205 lbs
Max. Speed	12 kmph / 7.5 mph
Weight Capacity	136 kg / 300 lbs
Ground Clearance	65 mm / 2.6"
Grade Climbable	10 degree
Curb Climbing	70 mm / 2.4"
Turning Radius	1400 mm / 55.1"
Suspension	Front & Rear
Brake	Electro-Mechanical
Seat Type	Sliding and Swivel Mid-Back
Seat Width	465 mm / 18.3"
Motor Size	700W, 4500 rpm
Battery Size	(2) 12V. 50Ah
Battery Weight	30.5 kg / 67.2 lbs
Travel Range	39 km / 24.2 miles (dependent on user weight and terrain)
Battery Charger	5A Off Board, Battery Level Indicator
Electronics	On /Off Key Switch, Speed Control Button
Colours Available	Blue, Burgundy



Your local dealer is: