

Karma Flexx Junior Wheelchair Script



Customer Details

Client							
Assessor							
Organisation							
	ACC			MOH			Other
User Weight							

Karma Flexx Junior Wheelchair Script

Please note: Scripts are not required for Karma Flexx Junior wheelchairs with the standard specifications (listed below). This script should only be used for clients with alternative requirements.

Standard Features:

- Standard seat and back upholstery
- 5" castor
- 20" pneumatic rear wheel (with quick release axle)
- Height adjustable, flip back armrest
- Heel loops
- 90 degree legrests with fixed angle foot rest
- Standard 2 point hip belt
- Back angle adjustment (-4° to 16°)
- Removable anti-tippers
- Push to lock brakes
- 80kg user weight limit

Wheelchair Type

<input type="checkbox"/>	Self Propelling
--------------------------	-----------------

- Quick Release 20" Pneumatic Spoked Rear Wheel
- Aluminium Handrim

<input type="checkbox"/>	Transit
--------------------------	---------

- Fixed 14" Solid Tyre Plastic Rim Rear Wheels

Seat Size Width x Depth

<input type="checkbox"/>	10 x 11
<input type="checkbox"/>	12 x 12
<input type="checkbox"/>	14 x 14

Front Seat Height

		Seat to Floor Height				
Castor Size	5"	15.5"	16"	17"		
	6"		16"	17"	18"	

Back Canes

<input type="checkbox"/>	Angle Adj.Std
<input type="checkbox"/>	Angle Adj Folding\$
<input type="checkbox"/>	RHS Height Adj Push Handle\$

Rear Seat Height

		Rear Seat Height						
Rear Wheel Size	14"	14.5"	15"					
	20"		15"	16"	17"	18"		
	22"			16"	17"	18"	19"	

Back Support

<input type="checkbox"/>	Standard (Tension adjustable with padding).....Std
<input type="checkbox"/>	Padded for Folding & Height Adj Push Canes.....\$
<input type="checkbox"/>	Omit Back Upholstery.....N/C
<input type="checkbox"/>	Alternative Back Support (please specify below).....\$

Alternative Back Support

Seat Upholstery

<input type="checkbox"/>	Standard (Tension adjustable with padding).....Std
<input type="checkbox"/>	Cushion (please specify below).....\$
<input type="checkbox"/>	Spare Cushion Cover.....\$

Rear Wheel Size

<input type="checkbox"/>	14" Solid Rim (Transit Option).....Std
<input type="checkbox"/>	20" Spoke.....Std
<input type="checkbox"/>	22" Spoke.....\$

Rear Wheel Option

<input type="checkbox"/>	Transit Tyre Plastic Rim (Available only as 14" Rear wheels).....Std
<input type="checkbox"/>	Pneumatic Tyre Spoked with Aluminium Handrim.....Std
<input type="checkbox"/>	Solid Tyre Spoked Rim with Aluminium Handrim.....\$

Leg Support

90° Swing-away	<input type="checkbox"/>	Left	<input type="checkbox"/>	Right	<input type="checkbox"/>	Std
<input type="checkbox"/>	Elevating Leg Rests.....\$					

Foot Support

Flip-up Angle & Depth Adjustable <input type="checkbox"/>	<input type="checkbox"/>	Left	<input type="checkbox"/>	Right	<input type="checkbox"/>	Std
---	--------------------------	------	--------------------------	-------	--------------------------	-----

Heel Loops

<input type="checkbox"/>	Heel LoopsStd
--------------------------	---------------------

Notes

Arm Support

<input type="checkbox"/>	Flip-back Height Adjustable.....Std
--------------------------	-------------------------------------

Wheel Locks

<input type="checkbox"/>	Push to Lock with Brake extension.....Std
--------------------------	---

Anti Tippers

<input type="checkbox"/>	Removable anti-tipper.....Std
--------------------------	-------------------------------

Accessories

<input type="checkbox"/>	Tray.....\$
--------------------------	-------------

Positioning Belts

<input type="checkbox"/>	Standard Hip BeltStd
<input type="checkbox"/>	Bodypoint 2 Point Belt S x 32cm Pads.....\$
<input type="checkbox"/>	Bodypoint 2 Point Belt S x 38cm Pads.....\$
<input type="checkbox"/>	Bodypoint 4 Point Belt S x 32cm Pads.....\$
<input type="checkbox"/>	Bodypoint 4 Point Belt S x 38cm Pads.....\$

Transit Option

<input type="checkbox"/>	Transit Tie Down Option.....\$
--------------------------	--------------------------------