

Account Number:

Date:

Address:

Quote: PO Number:

City:

Phone:

Associated Quote/Order#:

Email/Fax:

Effective: 5/22/2023

Catalyst 5 - Transit Included (K0005) \_\_\_\_\_ \$

Catalyst 5 - Transit Not Included (K0005) \_\_\_\_\_ \$

300 lb weight capacity

## Frame

### Seat Width

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	N/C
20" (E2201)	\$
21" (E2201)	\$
22" (E2201)	\$

### Seat Depth

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	\$75
20" (E2203)	\$

### Transit Option

Transit System (Weight Limit 300lbs)	\$
<i>Not available with Removable Seat Pan or 350 lb Weight Capacity</i>	
Omit	N/C

### Wheelchair Drive Ready

*Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Spinergy wheels or Height Adjustable Flip Back T-Arm.*

e-motion Ready	N/C
twion Ready	N/C

### Titanium Option

Titanium Frame \_\_\_\_\_ \$  
*Includes hangers and base frame. Does not include Cross Braces.*

### Weight Capacity

300 lbs	N/C
350 lbs (Heavy Duty)	\$
<i>Requires minimum 6" caster height</i>	

### Cross Brace

Standard Cross Braces	N/C
Seating System Cross Braces	N/C
<i>Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.</i>	
Titanium Cross Braces	\$
<i>Not available with 21" and 22" Seat Widths</i>	

### Seat Pan

Removable Aluminum Seat Pan E2231  
*1/2" Shorter than Depth* \_\_\_\_\_ \$

### Rear Frame Type

*No selection needed with Titanium Frame*

Standard (Curved)	N/C
<i>Select from Standard or Angle Adjustable Back Post on chart on the next page</i>	
Depth Adjustable (K0108)	\$
<i>Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.</i>	
Reclining Frame	\$
<i>Select from Reclining Back Post in chart on the next page. Not available with Seating Dynamics Rocker Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.</i>	

## Backrest

Depth Adjustable not available with Titanium Frame

Back Post	Price	Height Range	Back Height				
			8.5"	9.5"	10.5"	11.5"	
Standard	N/C	Short <sup>3,4</sup>	8.5"	9.5"	10.5"	11.5"	
		Medium <sup>4</sup>	12.5" <sup>6</sup>	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	
Angle Adjustable w/ Push Handle	\$	Short <sup>4</sup>	14"	15"	16"		
		Medium <sup>4</sup>	16"	17"	18"		
		Tall	18"	19"	20"		
Stroller Handle Back Post (Depth Adjustable)	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle (Depth Adjustable)	\$	Short <sup>1</sup>	13"	14"	15"	16"	17"
		Medium	16"	17"	18"	19"	20"
		Tall	20" <sup>8</sup>	21" <sup>5</sup>	22" <sup>5</sup>	23" <sup>5</sup>	24" <sup>5</sup>
Reclining w/ Adjustable Handle (E1225) <sup>1,7</sup>	N/C	Short	22"				
		Tall	26"				

1 Not available with Seating Dynamics Rocker Back

2 Not available with Standard Back Upholstery

3 Not available with Tension Adjustable Back Upholstery

4 Not available with Standard Back Upholstery in combination with 21" and 22" Seat Widths

5 Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths

6 Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width

7 Not available with O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole

8 Not available with Angle Adjustable Locking Extendable Flip Up Armrest

9 Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrests

### Reclining Back Options

Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.

- Standard Spring \_\_\_\_\_ N/C
- Heavy Duty Spring \_\_\_\_\_ N/C

### Standard Back Post

Select option below only if a Standard Back Post option is selected on chart above

- Straight with Push Handle \_\_\_\_\_ N/C
- Straight without Push Handle \_\_\_\_\_ N/C
- 8 Degree Bend with Push Handle \_\_\_\_\_ \$  
Not available with Short Back Height Range
- Straight with Fold Down Push Handle \_\_\_\_\_ \$  
Not available with 350 lb Weight Capacity

### Push Handle Options

Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 lb Weight Capacity

- Low Profile Push Handle \_\_\_\_\_ \$  
2" shorter than standard push handles

### Depth Adjustable Back Post Options

Options below may be selected only if a Depth Adjustable back post option is selected on the chart above

- Stroller Handle Extensions \_\_\_\_\_ \$  
Not available with Height Adjustable Straight with Push Handle Back Post
- Seating Dynamics Rocker Back (E2398) \_\_\_\_\_ \$  
Available only with T-Arms

### Back Upholstery Options

Not available with Reclining or Stroller Handle Back Post. If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Standard \_\_\_\_\_ N/C
- Tension Adjustable with Standard Cover (E2611) \_\_\_\_\_ \$
- Tension Adjustable with Air Mesh Breathable Cover (E2611) \_\_\_\_\_ \$
- Omit \_\_\_\_\_ (\$)

### Military Service Patches

Attached only if ordered with Tension Adjustable Back Upholstery

- Army Navy Marines Air Force Coast Guard US Purple Heart \_\_\_\_\_ \$
- Combo Army/Purple Heart Combo Navy/Purple Heart
- Combo Marines/Purple Heart Combo Air Force/Purple Heart
- Combo Coast Guard/Purple Heart \_\_\_\_\_ \$

## Frame Type / Footrest

### Frame Type

- Swing Away \_\_\_\_\_ N/C
- Super Low - Swing Away \_\_\_\_\_ N/C
- 70° Fixed Front \_\_\_\_\_ N/C
- 80° Fixed Front \_\_\_\_\_ N/C

### PRO ELR Covers

- Gel Pad Covers (Pair) K0108 \_\_\_\_\_ \$

### Footrest Hangers

Frame	Hanger	60°	70°	80°	90°	Retail
Swing Away	Extension Tube Mount	2		1		N/C
	Front Mount <sup>2</sup>			1	1	N/C
	Pro Elevating Leg Rests E0990					\$
	Omit Hangers					(\$)
Fixed	Extension Tube Mount			1		N/C

1 Not available with 7" or 8" casters

2 Not available with Titanium Frame

## Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length												
		11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"			
Extension Tube Hangers (60°,70°,80°)	Short (11" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"				
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				
Front Mount Hangers <sup>1</sup>	Standard Selections										Requires 3" Extension for \$35			
	Short Front Mount Hanger	4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"					10"	10.5"	11"
		7"	7.5"	8"	8.5"	9"	9.5"					11.5"	12"	12.5"
	Medium Front Mount Hanger	4" <sup>2</sup>	4.5" <sup>2</sup>	5"	5.5"	6"	6.5"	7"	7.5"			12"	12.5"	13"
		8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"			13.5"	14"	14.5"
	Long Front Mount Hanger	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"		14"	14.5"	15"
9.5"		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"		
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (12" - 17.5")	12 - 17.5												
	PRO ELR Low Pivot Setting Long (15" - 21")	15 - 21												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (10.5" - 15.5")	10.5 - 15.5												
	PRO ELR High Pivot Setting Long (13" - 19")	13 - 19												
Fixed Front Frame (70°,80°)	Short (11" - 15")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"					
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				

<sup>1</sup> Not available with Titanium Frame. See Minimum Front Seat Height chart below

<sup>2</sup> Available with 90° only

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

### 3" Extension

3" Extension	\$
3" Extension (Provided as Parts)	\$

### Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops (E0951)	\$
Adjustable Heel Loops (E0951)	\$
Omit Heel Loops	N/C

### Footplates

Composite	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable (K0040)	\$
<i>Not available with Front Mount Hanger or 14" Seat Width</i>	
Aluminum Angle Adjustable (K0040)	\$
Aluminum Locking Multi-Angle Adjustable Footplate (K0108)	\$
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable (K0040)	\$
<i>Not available with Front Mount Hanger or PRO ELR</i>	

### Hanger Release

If no selection is made Classic Hanger Release will be used

Classic	N/C
4-Way - In-Line	N/C
4-Way - Out	N/C

### Calf Straps

Velcro Adjustable (K0038)	\$
Padded Velcro Adjustable (K0038)	\$
Bodypoint Padded AeroMesh (K0038)	\$

### Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
				Smooth Cover Left	Air Mesh Cover Left	5x8 Pad Left Side	7x10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Left	Left Limb Support (E1020)	\$420	Omit Left Hanger	Smooth Cover Left	Air Mesh Cover Left	5x8 Pad Left Side	7x10 Pad Left Side	9x10 Pad Left Side	9x14 Pad Left Side	11x10 Pad Left Side	11x14 Pad Left Side
Right	Right Limb Support (E1020)	\$420	Omit Right Hanger	Smooth Cover Right	Air Mesh Cover Right	5x8 Pad Right Side	7x10 Pad Right Side	9x10 Pad Right Side	9x14 Pad Right Side	11x10 Pad Right Side	11x14 Pad Right Side

## Front Seat Height and Casters

### Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame			15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

### Casters

Caster Type		Price	Caster Size				
			4"	5"	6"	7"	8"
Caster Type	.75" Rollerblade	N/C					
	.75" Lighted Rollerblade	\$75					
	1" Poly	N/C					
	1" Poly Aluminum	\$115					
	1" Pneumatic	\$90					
	1.5" Poly	\$55					
	1.5" Soft Roll Aluminum	\$195					
	2" Poly	\$95					
	2" Pneumatic w/Foam Insert	\$240					

### Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

### Caster Housing

Standard Housing	N/C
Performance Housing	\$

Not available with 350 lb Weight Capacity

### Frog Legs Options

User Weight Under 80 lbs	N/C
User Weight 80-120 lbs	N/C
User Weight 120-160 lbs	N/C
User Weight 160-220 lbs	N/C
User Weight 220-260 lbs	N/C
User Weight 260-300 lbs	N/C

### Caster Pin Locks

Caster Pin Locks (K0073)	\$
--------------------------	----

Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

### Fork Options

Best Fit	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Standard Fork Stem	N/C
3/4" Longer Fork Stem	N/C
1 1/2" Longer Fork Stem	N/C
Frog Legs Suspension Fork	\$
<i>Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	
Frog Legs II Suspension System	\$
<i>Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	

## Wheel Locks

Wheel Lock	Wheel Lock Options
Push to Lock	N/C
Push to Lock Flush Mount	N/C
Pull to Lock	N/C
Low Profile Scissor Lock	\$
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>	
Drum Brake (E2228)	\$
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock (E2228)	\$
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock (E2228)	\$
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Hemi Wheel Lock - Left (K0108)	\$
<i>Extension Handle included</i>	
Hemi Wheel Lock - Right (K0108)	\$
<i>Extension Handle included</i>	
Omit	N/C
Extension Handles (E0961)	\$
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids (E0974)	\$
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

## Rear Seat Height and Wheels

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height									
		11.5" <sup>3</sup>	12" <sup>1</sup>	12.5" <sup>3</sup>	13" <sup>1</sup>	13.5" <sup>3</sup>	14"	14.5" <sup>3</sup>	15"	15.5" <sup>3</sup>	16" <sup>3</sup>
20"	Low	11.5" <sup>3</sup>	12" <sup>1</sup>	12.5" <sup>3</sup>	13" <sup>1</sup>	13.5" <sup>3</sup>	14"	14.5" <sup>3</sup>	15"	15.5" <sup>3</sup>	16" <sup>3</sup>
	Full	12" <sup>3</sup>	12.5" <sup>1</sup>	13" <sup>3</sup>	13.5" <sup>1</sup>	14" <sup>3</sup>	14.5"	15" <sup>3</sup>	15.5"	16" <sup>3</sup>	16.5" <sup>3</sup>
22"	Low	12.5" <sup>3</sup>	13" <sup>1</sup>	13.5" <sup>3</sup>	14" <sup>1</sup>	14.5" <sup>3</sup>	15"	15.5" <sup>3</sup>	16"	16.5" <sup>3</sup>	17" <sup>3</sup>
	Full	13" <sup>3</sup>	13.5" <sup>1</sup>	14" <sup>3</sup>	14.5" <sup>1</sup>	15" <sup>3</sup>	15.5"	16" <sup>3</sup>	16.5"	17" <sup>3</sup>	17.5" <sup>3</sup>
24"	Low	13.5" <sup>3</sup>	14" <sup>1</sup>	14.5" <sup>3</sup>	15" <sup>1</sup>	15.5" <sup>3</sup>	16"	16.5" <sup>3</sup>	17"	17.5" <sup>3</sup>	18" <sup>3</sup>
	Full	14" <sup>3</sup>	14.5" <sup>1</sup>	15" <sup>3</sup>	15.5" <sup>1</sup>	16" <sup>3</sup>	16.5"	17" <sup>3</sup>	17.5"	18" <sup>3</sup>	18" <sup>3</sup>
25"	Low	14" <sup>3</sup>	14.5" <sup>1</sup>	15" <sup>3</sup>	15.5" <sup>1</sup>	16" <sup>3</sup>	16.5"	17" <sup>3</sup>	17.5"	18" <sup>3</sup>	18.5" <sup>3</sup>
26"	Low		15" <sup>1,2</sup>	15.5" <sup>3</sup>	16" <sup>1</sup>	16.5" <sup>3</sup>	17"	17.5" <sup>3</sup>	18"	18.5" <sup>3</sup>	19" <sup>3</sup>
	Full	15" <sup>2,3</sup>	15.5" <sup>1,2</sup>	16" <sup>3</sup>	16.5" <sup>1</sup>	17" <sup>3</sup>	17.5"	18" <sup>3</sup>	18.5"	19" <sup>3</sup>	19.5" <sup>3</sup>

<sup>1</sup> Not available with Amputee Axle Plate

<sup>2</sup> Not available with Height Adjustable Flip Back T-Arm

<sup>3</sup> Requires Universal Axle Plate

One Arm Drive	Axle Plates
One Arm Drive - Quick Release (E0958)	No selection needed with One Arm Drive
Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.	Universal
Spoke wheel style is exclusive to OAD. Not available with 4 degrees of camber. Standard handrim will be Aluminum Anodized.	Amputee (E0959)
	Not available with Swing Away Armrests or Wheelchair Drive Ready Option.
	Required with Reclining Backrest.
Drive Side	
Left Side Drive	N/C
Right Side Drive	N/C
Drive Handrim	
Aluminum Drive Handrim	N/C
Poly Coated Drive Handrim	\$

## Rear Wheel Type

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Maxx Spoke	\$210					
	Maxx Spoke	N/C					
	Mag	N/C					
	Superlight Spoke <sup>1</sup>	\$275				<sup>2</sup>	
	Spinergy Spox	\$915					
	Spinergy LX	\$1148					

<sup>1</sup> Not available with 350 lb Weight Capacity

<sup>2</sup> Available only with Superlight and Plastic Coated Handrims

## Spinergy Options

Select spoke color and handrim tabs for Spinergy wheels

Black Spokes   Red Spokes   Blue Spokes   White Spokes \_\_\_\_\_ N/C

## Handrims

No selection needed with One Arm Drive

### Standard Handrims

Aluminum Anodized	N/C
Superlight Aluminum Anodized	\$
<i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels.</i>	
Plastic Coated	\$
Projection (8 Vertical Projections)	\$
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

### Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Lite Regular Grip	\$
Standard Regular Grip	\$
Lite Extra Grip	\$
Standard Extra Grip	\$
Lite Omit Grip	\$
Standard Omit Grip	\$

## Handrim Options

### Handrim Mount

If no selection is made, Standard Mount will be used.

Close Mount	\$30
-------------	------

### 3M Micro Grip

3M Micro Grip	\$
<i>Not available on Plastic Coated or Projection Handrims</i>	

## Rear Tire Type

### Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" (E2218)	N/C
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$
<i>Not available on Superlight wheel or 25" wheel size</i>	
Schwalbe HS-127 1 3/8" (E2211, E2212)	\$
<i>Available only on 24" wheel size. Not available on Mag wheel.</i>	
Schwalbe Downtown 1 3/8" (E2211, E2212)	\$

### Low Profile

Low Poly 1" (E2218)	N/C
<i>Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.</i>	
Sentinel Solid Tire 1" (E2218)	\$
<i>Not available on Superlight wheel or 20" and 22" Mag wheels</i>	
High Pressure 1" (E2211, E2212)	\$
<i>Not available on Mag wheel</i>	
Schwalbe One 1" (E2211, E2212)	\$
<i>Not available on Mag wheel, 20", 22" or 26" wheel sizes</i>	
Schwalbe Marathon Plus Evolution 1" (E2211, E2212)	\$
<i>Not available on Mag wheel, 20" or 22" wheel sizes</i>	

## Axles

No selection needed with One Arm Drive or Drum Brake

Quick Release	N/C
Quad Quick Release	\$
<i>Not available with Mag or Spinergy wheels</i>	

## Axle Position

### Rear Wheel Camber

0°   2°   4°	N/C
--------------	-----

### Center of Gravity

Rear   Center   Forward	N/C
-------------------------	-----

No selection needed with Amputee Axle Plate

## Armrests

Armrest Type	Length	
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$
	Full	
Height Adjustable T-Arm Low (6.5"-10.5") E0973 <sup>1</sup>	Desk	\$
	Full	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") E0973 <sup>2</sup>	Desk	\$
	Full	
Height Adjustable Flip Back T-Arm (9.5"-13.5") E0973 <sup>2</sup>	Desk	\$
	Full	

<sup>1</sup> If requested configuration is unavailable, Standard Height T-Arm will be used

<sup>2</sup> Available only with Curved Rear Frame

<sup>3</sup> Not available with Seat Depths less than 16"

<sup>4</sup> Available only with Standard Back Posts

<sup>5</sup> Available only with Depth Adjustable and Reclining Back Posts

<sup>6</sup> Not available with Arm Pads

<sup>7</sup> Not Available with 15" Rear Seat Height and 24" wheels

Armrest Type	Length	Price
Flip back (10.5") <sup>3,4,7</sup>	Desk	N/C
	Full	
Angle Adjustable Locking Extendable Flip Up E0973 <sup>5</sup>	Desk	\$
	Full	
Swing Away (10.5"-12.5") E0973 <sup>4,6</sup>		\$
Swing Away with Side Guard (10.5"-12.5") E0973 <sup>4,6</sup>		\$
Omit Armrests		N/C

### Arm Pads

Standard Arm Pad	N/C
Waterfall Arm Pad	\$
<i>Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low</i>	
Foam Grip	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

## Accessories

### Removable Side Guard

*Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests*

Adult Composite Side Guards (K0108)	\$
Pediatric Composite Side Guards (K0108)	\$
Aluminum Straight - Black Matte (K0108)	\$
Aluminum Straight - Match Frame Color (K0108)	\$
<i>Not available with Hydrographics or Titanium Frame</i>	
Aluminum Adjustable Fender w/Foam Pad - Black Matte (K0108)	\$
Aluminum Adjustable Fender w/Foam Pad - Match Frame Color (K0108)	\$
<i>Not available with Hydrographics or Titanium Frame</i>	
Carbon Fiber Adjustable Fender (K0108)	\$

### Anti-tips

Rear Anti-Tips (E0971)	\$
Swing-In Anti-Tip Single - Left (E0971)	\$
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.</i>	
Swing-In Anti-Tip Single - Right (E0971)	\$
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.</i>	
Swing-In Anti-Tip Pair (E0971)	\$
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with Cane and Crutch Holder or O2 Holder.</i>	
Omit	N/C

### Accessories

#### Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
2" Airline Buckle Pelvic Positioning Belt (E0978)	\$
1" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$
Small Bodypoint Evoflex (16-24") (E0978)	\$
Medium Bodypoint Evoflex (19-29") (E0978)	\$
Large Bodypoint Evoflex (25-37") (E0978)	\$
Belt Mounting Band Clamp	\$

#### O2 Holder

*Not Available with Depth Adjustable Rear Frame*

O2 Holder - Left Side (E2208)	\$
O2 Holder - Right Side (E2208)	\$

#### Additional Accessories

Tool Kit	\$
Removable Underseat Pouch	\$
Pediatric Ki Mobility Backpack	\$
Ki Mobility Backpack	\$
Luggage Carrier	\$
Black Neoprene Impact Guards	\$
Clear Spoke Guards (K0065)	\$
Black Spoke Guards (K0065)	\$
Cane and Crutch Holder	\$
IV Pole	\$

## Color

### Gloss Finish

Black	N/C
Blue Skies	N/C
Bubblegum Pink	N/C
Candy Blue	N/C
Candy Purple	N/C
Candy Red	N/C
Copper	N/C
Granny Smith	N/C
Metallic Orange	N/C
Yellow	N/C

### Matte Finish

Military Green	N/C
Charcoal Gray	N/C
Emerald Green	N/C
Iridescent Burgundy	N/C
Iridescent Navy	N/C
Retro Red	N/C
Royal Blue	N/C
Shadow Black	N/C
Smooth White	N/C

### Hydrographics

AOR1 (Desert Digital Camo)	\$
NEXT Camo	\$
Pink Camo	\$
Carbon Fiber	\$
<i>Select Hydrographics Base Color</i>	
Muddy Girl	\$
Zebra	\$
<i>Select Hydrographic Base Color. Zebra not available in Black.</i>	

### Hydrographics Base Color

*Required with Hydrographics*

Black	Blue	Pink	Granny Smith	N/C
Purple	Red			

### Colored Anodizing Package

*Select desired color. If no selection is made, black will be used.*

	Option	Standard Caster Housing <sup>2</sup>	Caster Fork <sup>3</sup>	Axle Plate <sup>4</sup>
		\$	\$	\$
Color	Black <sup>1</sup>			
	Blue			
	Bright Green			
	Orange			
	Pewter			
	Pink			
	Purple			
	Red			

<sup>1</sup> Black is no charge

<sup>2</sup> Available only if corresponding option selection was made

<sup>3</sup> Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System

<sup>4</sup> Color Anodized not available with amputee axle plate



Catalyst Front Seat Height Matrix																			
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster			
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4	
15.5 <sup>1</sup>	14 <sup>1</sup>			17.5	16	1 1/2		18.5	17	1 1/2		18.5	17.5	3/4		20.5	n/a	1 1/2	
16	14.5	3/4	4"	16	14.5	Std	5"	17	15.5	Std	6"	19	18	1 1/2	7"	19	17.5	Std	7"
16.5 <sup>1</sup>	15 <sup>1</sup>			16.5	15			17.5	16			18	16.5			19.5	18		
17	15.5	1 1/2	4"	16.5	15.5	3/4	5"	18	16.5	Std	6"	20	18.5	Std	7"	20	18.5	Std	7"
15.5	14			17	16			17.5	16			18	17			19	17.5		
16	14.5	Std	5"	17.5	16	1 1/2	6"	18	17	3/4	7"	19	17.5	Std	8"	19.5	18.5	3/4	8"
16	15			18	16.5			18.5	17.5			19.5	18			20	18.5		
16.5	15.5	3/4	5"	16.5	15	Std	6"	18.5	17	1 1/2	7"	20	18.5	3/4	8"	20.5	n/a	1 1/2	8"
17	15.5			17	15.5			19	17.5			19.5	18			21	n/a		
17.5	16	1 1/2	5"	17.5	16	Std	6"	19.5	18	1 1/2	7"	20.5	19	Std	8"	21	n/a	1 1/2	8"
16.5	15			17.5	16			18.5	17			19.5	18			20.5	19		
16.5	15	Std	6"	17.5	16	3/4	7"	18	16.5	Std	8"	20	19	1 1/2	9"	20.5	n/a	1 1/2	9"
17	15.5			18	16.5			18.5	17			19.5	18			20.5	n/a		
17	16	3/4	6"	18.5	17	Std	7"	19	17.5	Std	8"	20	18.5	1 1/2	9"	20.5	19	1 1/2	9"
17.5	16.5			18	16.5			19.5	18			20.5	19			21	n/a		
18	16.5	1 1/2	6"	18.5	17	1 1/2	7"	18.5	17.5	3/4	8"	21	n/a	Std	9"	21	n/a	1 1/2	9"
18.5	17			19	17.5			19.5	18			20.5	19			21	n/a		
18.5	17.5	3/4	7"	19	17.5	Std	8"	19.5	18	1 1/2	9"	20.5	19	Std	10"	21	n/a	1 1/2	10"
19	18			18.5	17.5			20	18.5			20.5	19			21	n/a		
19.5	18.5	1 1/2	7"	19	18	3/4	8"	20.5	19	1 1/2	9"	21	n/a	Std	10"	21	n/a	1 1/2	10"
20	19			19.5	18.5			20.5	19			21	n/a						
				19.5	18	1 1/2	8"	6" x 2" caster not available with 5" fork											
				20	18.5														
				20.5	19														

<sup>1</sup> Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

## Cushions and Backrests

**AXIOM**

**Ki Mobility**

**ADI**  
A Stealth Product

**VARILITE**

**R I O E**

Cushion:

Backrest:

Other Accessories: