



Company: _____

Therapist: _____ Date: _____

☐ Acc ☐ MoH

☐ Quote ☐ PO Number: _____

Marked For: _____

Rogue K0005, E1235 \$
Rogue Standard Frame. 125 kg weight capacity.

Ship To: _____

Address: _____

City: _____ State: _____

Phone: _____

Email/Fax: _____

Notes: _____ Effective 23/04/2019

Frame

Seat Width / Seat Depth

Select box corresponding to desired Seat Width and Depth.
See references A and B on page 7.

		Seat Width								
		12"	13"	14"	15"	16"	17"	18"	19"	20"
		Price HCPCS	N/C	N/C	N/C	N/C	N/C	N/C	N/C	\$ E2201
Seat Depth	14"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	15"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	17"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	18"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19"	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	\$ E2203	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Transit Option (Tie Downs)

Cannot be retrofitted

☐ Transit System K0108 \$

Backrest

Back Post / Height Range / Back Height

Select Back Height in row corresponding to desired Back Post and Height Range.

Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	Short	<input type="checkbox"/> 11 ²	<input type="checkbox"/> 12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	
		Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		Tall	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	
Fixed Height ¹	N/C		<input type="checkbox"/> 9 ²	<input type="checkbox"/> 10 ²	<input type="checkbox"/> 11 ²	<input type="checkbox"/> 12	<input type="checkbox"/> 13
			<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18

¹ Not available with back upholstery

² Not available with Fold Down Push Handle

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 7.

☐ 0" (Performance) N/C
☐ 1" N/C
☐ 2" (Classic) N/C

Adjustable Depth

Adds the capability to grow seat in depth

☐ No Rear Adjustment N/C
☐ 1/2" N/C
☐ 1" \$
☐ 2" \$

Seat Pan

☐ Solid Seat Pan E2231 \$

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 7.

☐ 90° \$
☐ 85° N/C
☐ 80° N/C
☐ 75° N/C
☐ 70° N/C

Push Handle

☐ Omit Push Handle N/C
☐ Standard Push Handle \$
☐ Low Profile Push Handle \$
 2" shorter than standard push handles
☐ Bolt-On Push Handle \$
 Handles will be 1" - 2" above back height
☐ Fold Down Push Handle \$

Stroller Handle

Stroller Handles are removable and height adjustable

☐ Single Stroller Handle \$
☐ Double Stroller Handle \$
 Not available with Push Handles

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0025.30 Pricing and specifications subject to change without notice

Page 1 of 7

Backrest continued

Rigidizer Bar

☐ Adjustable Height Rigidizer Bar _____ \$

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

☐ Standard ☐ +3° ☐ +6° ☐ +9°
☐ -3° ☐ -6° ☐ -9° ☐ -12°

Back Upholstery

Available only with Height Adjustable Back Posts

Padded Velcro Adjustable

☐ Standard Cover _____ Std

☐ Air Mesh Breathable Cover **New!** _____ \$

Tension Adjustable E2611

☐ Standard Cover _____ \$

☐ Air Mesh Breathable Cover **New!** _____ \$

☐ Omit Back Upholstery _____ N/C

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height in row corresponding to Footrest Taper and Caster Size. 6" casters not available with 85° and 90° Front Frame Bend.

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube.

2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 7.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 7.

Footrest Taper	Caster Size	Front Seat Height														
Straight	3"	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5
	4"	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
	5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5		
	6"	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5			
1" Taper	3"	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5						
	4"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							
	5"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5								
	6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5									
2" Taper	3"	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5						
	4"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5							
	5"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5								
	6"	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5									

Caster Type

☐ .75" Rollerblade _____ N/C

Available only with 3" and 4" Caster Sizes

☐ .75" Lighted Rollerblade _____ \$

Available only with 3" and 4" Caster Sizes

☐ 1" Poly _____ N/C

Not available with 3" Caster Size

☐ 1" Poly Aluminum _____ \$

Available only with 4" and 5" Caster Sizes

☐ 1" Pneumatic _____ \$

Available only with 6" Caster Size

☐ 1.5" Poly _____ \$

Not available with 3" Caster Size

☐ 1.5" Soft Roll Aluminum _____ \$

Fork Options

☐ Standard Fork _____ Std

☐ Single Sided Fork _____ \$

☐ Frog Legs Ultra-Sport Fork - Black E1015 _____ \$

Not included in Colored Anodizing Package

End User Weight Required: _____ kg

☐ Frog Legs PHASE 2 Carbon Fork _____ \$

Lightweight shock absorbing fork. Not included in Colored Anodizing Package.

Footrest

Footplate

- ☐ Tubular "Open" Footplate _____ Std
- ☐ Tubular Footplate with ABS Cover _____ \$30
Not available w/footrest widths greater than 15.5". See Reference G on page 7.
- ☐ Aluminum Angle Adjustable K0040 _____ \$95
- ☐ Aluminum Angle Adjustable with Riser K0040, K0037 _____ \$165
- ☐ High Mount Aluminum Angle Adjustable K0040, K0037 _____ \$165
- ☐ Angle Adjustable Flip Under K0040, K0037 _____ \$185
Not available with footrest widths less than 9.5". See Reference G on page 7.
- ☐ High Mount Angle Adjustable Flip Under K0040, K0037 _____ \$270
Not available w/footrest widths greater than 15.5". See Reference G on page 7.

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest.
Footrest length at least 2.5" shorter than Front Seat Height
recommended. See reference H on page 7.

<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"
<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"
<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"
<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"

Footrest Protector

- ☐ Footrest Protector _____ N/C
Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels.
Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 7.

Wheel Size	Rear Seat Height															
20"	<input type="checkbox"/> 12.5 ¹	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
22"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20		
24"	<input type="checkbox"/> 14.5 ¹	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20				
25"	<input type="checkbox"/> 15 ¹	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20					
26"	<input type="checkbox"/> 15.5 ¹	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20						

¹ Available only with Low Profile tires

One Arm Drive

- One Arm Drive E0958 _____ \$
Available only with 24" Mag wheel. Axle not quick release. Standard Handrim will be Aluminum Anodized.
- Drive Handrim Side: ☐ Left ☐ Right
- Drive Handrim Type: ☐ Aluminum _____ N/C
☐ Poly Coated _____ \$

Rear Wheel Type

- ☐ MAXX Performance Spoke _____ N/C
- ☐ Mag _____ N/C
Not available in 25" or 26" wheel sizes
- ☐ Superlight Spoke _____ \$
Available only in 24" and 25" wheel sizes. 25" wheel size available only with Superlight and Plastic Coated Handrims.

- Spinergy Spox _____ \$
Not available in 20" wheel size
- Color: ☐ Black ☐ Red ☐ White ☐ Blue ☐ Yellow

- Handrim: ☐ Short Tab ☐ Long Tab
Available only with Aluminum Anodized and Plastic Coated Handrims

- Spinergy LX _____ \$
Not available in 20" wheel size

- Color: ☐ Black ☐ Red ☐ White ☐ Blue ☐ Yellow
☐ Pink ☐ Green ☐ Orange ☐ Purple **New!**

- Handrim: ☐ Short Tab ☐ Long Tab
Available only with Aluminum Anodized and Plastic Coated Handrims

Rear Tire Type

Full Profile

- ☐ Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
- ☐ Full Poly 1 3/8" E2218 _____ \$
Not available on Superlight wheel or 25" or 26" wheel sizes
- ☐ Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$
Not available on Superlight wheel or 25" wheel size
- ☐ Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$
Available only on 24" wheel size. Not available on Mag wheel.
- ☐ Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$
Available only on 24" wheel size
- ☐ Mountain 2" Black E2211, E2212 _____ \$
Available only on 22", 24" and 25" Spoke and Spinergy wheels

Low Profile

- ☐ Sentinel Solid Tire 1" E2218 _____ \$
Not available on Superlight wheel or 20" and 22" Mag wheels
- ☐ High Pressure 1" E2211, E2212 _____ \$
Not available on Mag wheel
- ☐ Schwalbe One 1" E2211, E2212 _____ \$
Available only on 24" and 25" wheel sizes. Not available on Mag wheel.
- ☐ Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$
Not available on Mag wheel, 20" or 22" wheel sizes

Rear Seat Height and Wheels continued

Handrim

No selection needed with One Arm Drive

Standard Handrims

- ☐ Aluminum Anodized _____ Std
- ☐ Superlight Aluminum Anodized _____ \$
Available only on 24" and 25" wheel sizes. Available only on MAXX Performance, Superlight or Spinergy wheels.
- ☐ Plastic Coated _____ \$
- ☐ Projection (8 Vertical Projections) _____ \$
Available only with full profile tires in combination with MAXX Performance, 22"-26" Spoke, or 24" Mag wheel
- ☐ Spinergy Flex Rim _____ \$
Available only on 24" and 25" Spinergy LX Wheels
- ☐ Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- ☐ Natural Fit LT Regular Grip _____ \$
- ☐ Natural Fit Standard Regular Grip _____ \$
- ☐ Natural Fit LT Super Grip _____ \$
- ☐ Natural Fit Standard Super Grip _____ \$
- ☐ Natural Fit LT Omit Grip _____ \$
- ☐ Natural Fit Standard Omit Grip _____ \$

Handrim Option

- ☐ 3M Micro Grip _____ \$
Not available on Plastic Coated or Projection Handrims

Axle

- ☐ Quick Release _____ Std
- ☐ Quad Quick Release _____ \$
Not available with Mag or Spinergy wheels

Wheel Locks

Wheel Lock

- ☐ Push to Lock _____ Std
- ☐ Push to Lock Flush Mount _____ N/C
- ☐ Pull to Lock _____ N/C
- ☐ Under Seat Scissor Lock _____ \$
- ☐ Low Profile Scissor Lock _____ \$
Not available with Low Poly, Full Poly or Sentinel Solid tires
- ☐
- ☐ Omit Wheel Locks _____ N/C

Rear Wheel Camber

No selection needed with One Arm Drive

- ☐ 0° ☐ 2° ☐ 4° ☐ 6° ☐ 8°

Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable Armrests. 1" minimum Rear Wheel Spacing with 0° Rear Wheel Camber. See reference K on page 7. No selection needed with One Arm Drive or Wheelchair Drive Ready options.

- ☐ .75" ☐ 1" ☐ 1.25" ☐ 1.5" ☐ 1.75"

Center of Gravity Preset

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 7.

Amputee Positions

<input type="checkbox"/> -1.5"	<input type="checkbox"/> -1.25"	<input type="checkbox"/> -1"	<input type="checkbox"/> -.75"
<input type="checkbox"/> -.5"	<input type="checkbox"/> -.25"	<input type="checkbox"/> 0"	<input type="checkbox"/> .25"

Standard Positions

<input type="checkbox"/> .5"	<input type="checkbox"/> .75"	<input type="checkbox"/> 1"	<input type="checkbox"/> 1.25"
<input type="checkbox"/> 1.5"	<input type="checkbox"/> 1.75"	<input type="checkbox"/> 2"	<input type="checkbox"/> 2.25"
<input type="checkbox"/> 2.5"	<input type="checkbox"/> 2.75"	<input type="checkbox"/> 3"	<input type="checkbox"/> 3.25"
<input type="checkbox"/> 3.5"	<input type="checkbox"/> 3.75"	<input type="checkbox"/> 4"	<input type="checkbox"/> 4.25"

Wheelchair Drive Ready

Ki Chairs will fit most power assist options

Wheel Lock Options

- ☐ Extension Handles E0961 _____ \$
Available only with Push to Lock or Pull to Lock Wheel Locks
- ☐ Grade Aids E0974 _____ \$
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Armrests

Armrest		Please select Arm Pad option	
<input type="checkbox"/>	Height Adjustable T-Arm (8.5" - 13.5") E0973	\$	
<input type="checkbox"/>	Height Adjustable T-Arm Low (7.5" - 11.5") E0973	\$	
Arm Pad			
Not available with Swing Away Armrests			
<input type="checkbox"/>	Desk	<input type="checkbox"/>	Full

Armrest continued		
<input type="checkbox"/>	Swing Away (8.5" - 12.5") E0973 <i>Not available with Arm Pads</i>	\$
<input type="checkbox"/>	Swing Away with Side Guard (8.5" - 12.5") E0973 <i>Includes direct mount receiver. Not available with Arm Pads.</i>	\$
<input type="checkbox"/>	Omit Armrests	N/C

Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Direct Mount		Clamp-On Mount	
<input type="checkbox"/> Composite Side Guards K0108	\$	<input type="checkbox"/> Composite Side Guards K0108	\$
Aluminum Straight K0108		Aluminum Adjustable Fender with Foam Pad K0108	
<input type="checkbox"/> Black Matte	\$	<input type="checkbox"/> Black Matte	\$
<input type="checkbox"/> Match Frame Color	\$	<input type="checkbox"/> Match Frame Color	\$
<i>Not available with Hydrographics</i>		<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender with Foam Pad K0108		<input type="checkbox"/> Carbon Fiber Adjustable Fender K0108	\$
<i>Not available with Adjustable Depth greater than 1/2"</i>			
<input type="checkbox"/> Black Matte	\$		
<input type="checkbox"/> Match Frame Color	\$		
<i>Not available with Hydrographics</i>			
<input type="checkbox"/> Carbon Fiber Adjustable Fender K0108	\$		
<i>Not available with Adjustable Depth greater than 1/2"</i>			

Accessories

Anti-Tip		Belt Mounting Kit	
<input type="checkbox"/> Rear Anti-Tips E0971	\$	<input type="checkbox"/> Ki Mobility Belt Mounting Clamp	\$
		<input type="checkbox"/> Bodypoint Belt Mounting Bracket Kit	\$
Calf Strap		Other Accessories	
<input type="checkbox"/> Velcro Adjustable K0038	\$	<input type="checkbox"/> Tool Kit	\$
<input type="checkbox"/> Padded Velcro Adjustable K0038	\$	<input type="checkbox"/> Removable Underseat Pouch	\$
<input type="checkbox"/> Bodypoint Padded AeroMesh K0038	\$	<input type="checkbox"/> Ki Mobility Backpack	\$
		<input type="checkbox"/> Luggage Carrier	\$
Positioning Belt		<input type="checkbox"/> Black Neoprene Impact Guards	\$
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$	<input type="checkbox"/> Under Chair Carrier and Removable Bag	\$
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$	<i>Not available with Seat Widths less than 14"</i>	
<input type="checkbox"/> 1" Auto Buckle E0978	\$	<input type="checkbox"/> Clear Spoke Guards K0065	\$
<input type="checkbox"/> 1" Padded Auto Buckle E0978	\$	<i>Not available with Mag wheel</i>	
<input type="checkbox"/> 2" Airline Buckle E0978	\$	<input type="checkbox"/> Black Spoke Guards K0065	\$
		<i>Not available with 20" wheel size or Mag wheel</i>	
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$	<input type="checkbox"/> Cane and Crutch Holder	\$
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")			

Color

Gloss Frame Finish

<input type="checkbox"/> Black	N/C
<input type="checkbox"/> Blue Skies	N/C
<input type="checkbox"/> Bubblegum Pink	N/C
<input type="checkbox"/> Candy Blue	N/C
<input type="checkbox"/> Candy Purple	N/C
<input type="checkbox"/> Candy Red	N/C
<input type="checkbox"/> Copper	N/C
<input type="checkbox"/> Granny Smith	N/C
<input type="checkbox"/> Metallic Fuchsia New!	N/C
<input type="checkbox"/> Metallic Orange	N/C
<input type="checkbox"/> Silver	N/C
<input type="checkbox"/> Yellow	N/C

Matte Frame Finish

<input type="checkbox"/> Military Green	N/C
<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Emerald Green	N/C
<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Retro Red	N/C
<input type="checkbox"/> Royal Blue	N/C
<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Smooth White	N/C

Hydrographics

Visit www.kimobility.com for current Hydrographic options. Note selection below

☐ _____ \$
Not available with Transit Option

Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

	Option	Back Plates ²	Camber Tube ³	Towers	CG Clamps	Side Guard Receivers	RP Tuning Housing	Front Caster Forks ⁴
	Retail	\$	\$	\$	\$	\$	\$	\$
Color	Black ¹	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Bright Green	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Gold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Orange	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pewter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pink	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Purple	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹ Black is no charge

² Not available with Transit Option

³ Not available with Wheelchair Drive Ready

⁴ Available only with Standard or Single Sided

Cushions and Backrests

AXIOM

ottobock.

KiMobility

ADI
A Stealth Product

VARILITE

R I O E

Cushion:

Backrest:

References

A. Seat Width *See Diagram 1*

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth *See Diagram 3*

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position *See Diagram 2*

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1° or 2° seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend *See Diagram 3*

Angle between front seat tube and ground.

E. Footrest Taper *See Diagram 1*

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height *See Diagram 1*

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width *See Diagram 1*

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
12"	9.5"	7.5"	
13"	10.5"	8.5"	
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"

H. Seat to Footrest Length *See Diagram 3*

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height *See Diagram 3*

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset *See Diagram 3*

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

K. Rear Wheel Spacing *See Diagram 1*

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2




Diagram 3

