

karma

S-ERGO

The Karma S-Ergo Wheelchair is a very lightweight self-propelling wheelchair, designed to fit the curves of your body. It folds easily, is very lightweight and extremely compact, making it very easy for clients and caregivers to manage.



**ALLIED
MEDICAL**
LIMITED



Performance, Style and Comfort

+ UNIQUE ERGONOMIC SEATING

The S-Ergo seating system features a patented ergonomic design which acts to relieve pressure as well as promote good posture. The S-shaped seat fits with the natural curves of the body and works to reduce the user sliding forward. The extremely light frame makes it easy to lift in and out of vehicles, reducing any additional strain on carers.



+ FEATURES



- Flip back armrests
- Extended armrests that allow for wider seat area while minimising the overall width of the wheelchair
- Ergonomic shaped armrest designed to offer comfort and support
- Ergo push rim that aids wheelchair propulsion
- Swing-away footrests that swing both inwards and outwards
- Quick release, solid flat-free rear wheels
- Foldable back canes
- Attendant brakes
- Easily removable cushion from upholstery for cleaning purposes

+ ACCESSORIES

- Spokeguards
- Heel loops
- Stump support
- Folding headrest

+ SPECIFICATIONS



Seat Width
16" (40.5 cm)
18" (45.5 cm)




Seat Depth
17" (43.2 cm)




Wheels
7" solid (front)
24" solid (rear)



Seat to Floor Height
20"



Overall Width
16" Width: 590 mm
18" Width: 640 mm
20" Width: 690 mm



Folded Width
325 mm



Back Height
410 mm



Overall Length
14" Wheels: 910 mm
24" Wheels: 1020 mm



Overall Height
905 mm



Total Weight
Approx 13.2 kg



User Weight Limit
115 kg



**ALLIED
MEDICAL
LIMITED**

0800 31 61 81
www.alliedmedical.co.nz