

***Weights  
just  
8.5 kg!***

**KARMA**

Karma KM2501 Ergolite

**Comfortable and light weight, with all  
the features you expect from a folding  
transit wheelchair, and more!**



***Ultra-light  
Folding  
Transit  
Wheelchair***



**ALLIED  
MEDICAL**  
LIMITED

## MOBILITY HAS NEVER BEEN SO EASY

This amazingly light and compact transit wheelchair is ideal for outings and travellers. It folds down to take up virtually no space in the boot of a car and weighs only 8.5 kg – making it easy for anyone to lift into a vehicle.

The patented S-Ergo seat system is designed to fit the natural curves of the human body that creates a comfortable seating position that promotes pressure redistribution and reduces the possibility of sliding forward.

Available in sizes of 16" wide x 16" deep, and 18" wide x 16" deep, the KM2501 is the perfect chair if you're seeking a solution for outings to the park, shopping mall, or anywhere in the community for that matter!

It also has an optional carry bag that makes it easy to travel with – the world is your oyster!



## FEATURES

- Ergonomically designed armrests that allow you to get right up to a table or desk
- Calf strap
- Dual brake system – attendant brake as well as manual brakes on each wheel
- Folding backrest
- Pump and tool kit for ongoing standard maintenance
- **aegis** MICROBE SHIELD™ Microbe shield upholstery protects against odour, stains, bacteria, fungus and other microorganisms

## ACCESSORIES

- Quick release rear wheels
- Single wheeled antitipper
- Travel bag
- Positioning hip belt



Keep your Ergolite safe with a travel bag

## SPECIFICATIONS



**Seat widths**  
16" (40 cm),  
18" (45 cm)




**Seat depth**  
16.5" (42 cm)



**Wheels**  
6" (front)  
14" (rear)



**Seat to floor height**  
19"



**Overall width**  
57 cm (16" w),  
62 cm (18" w)



**Folded width**  
27 cm



**Back height**  
16.5" (42 cm)



**Overall length**  
94 cm



**Overall height**  
91 cm



**Total weight**  
8.5 kg (16" w)  
8.5 kg (18" w)



**User weight limit**  
100 kg



**0800 31 61 81**  
[www.alliedmedical.co.nz](http://www.alliedmedical.co.nz)