

Q-LOGIC

DRIVE CONTROL SYSTEM

SEAT CONFIGURATION SETTINGS (AAM)

CONFIGURATION NUMBER	A1 TILT	A2 RECLINE	A3 ELEVATE & LIFT	A4 LEFT LEG	A4 BOTH LEGS	A4 FOOTPLATFORM	A5 RIGHT LEG	A5 BACK	FUNCTIONS ACTIVATED BY CONFIGURATION
1				X			X		INDEPENDENT LEGS
2					X				COMBINED LEGS
3						X			ARTICULATING FOOT PLATFORM
4			X						ELEVATING SEAT
5			X	X			X		ELEVATING SEAT W/ INDEPENDENT LEGS
6			X		X				ELEVATING SEAT W/ COMBINED LEGS
7			X			X			ELEVATING SEAT & ARTICULATING FOOT PLATFORM
8		X							RECLINE
9		X		X			X		RECLINE/INDEPENDENT LEGS W/ SYNCHRONIZED RECLINE & LEGS
10		X			X				RECLINE /COMBINED LEGS W/ SYNCHRONIZED RECLINE & LEGS
11		X				X			RECLINE/ ARTICULATING FOOT PLATFORM W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
12		X			X			X	RECLINE/ COMBINED LEGS/ BACK W/ SYNCHRONIZED RECLINE & LEGS
13		X				X		X	RECLINE/ ARTICULATING FOOTPLATFORM/ BACK W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
14		X	X						RECLINE & ELEVATING SEAT
15		X	X	X			X		RECLINE/ELEVATING SEAT/ INDEPENDENT LEGS W/ SYNCHRONIZED RECLINE & LEGS
16		X	X		X				RECLINE/ELEVATING SEAT/ COMBINED LEGS W/ SYNCHRONIZED RECLINE & LEGS
17		X	X			X			RECLINE/ELEVATING SEAT/ FOOT PLATFORM W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
18		X	X		X			X	RECLINE/ELEVATING SEAT COMBINED LEGS/BACK W/ SYNCHRONIZED RECLINE & LEGS
19		X	X			X		X	RECLINE, ELEVATING SEAT, FOOT PLATFORM, BACK W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
20	X								TILT
21	X			X			X		TILT & INDEPENDENT LEGS
22	X				X				TILT & COMBINED LEGS
23	X					X			TILT & ARTICULATING FOOT PLATFORM
24	X		X						TILT & ELEVATING SEAT
25	X		X	X			X		TILT/ELEVATING SEAT/& INDEPENDENT LEGS
26	X		X		X				TILT/ELEVATING SEAT/ & COMBINED LEGS

SEAT CONFIGURATION SETTINGS (AAM)

CONFIGURATION NUMBER	A1 TILT	A2 RECLINE	A3 ELEVATE & LIFT	A4 LEFT LEG	A4 BOTH LEGS	A4 FOOTPLATFORM M	A5 RIGHT LEG	A5 BACK	FUNCTIONS ACTIVATED BY CONFIGURATION
27	X		X			X			TILT/ELEVATING SEAT/& ARTICULATING FOOT PLATFORM
28	X	X							TILT & RECLINE
29	X	X		X			X		TILT/RECLINE/INDEPENDENT LEGS W/ SYNCHRONIZED RECLINE & LEGS
30	X	X			X				TILT/RECLINE/ COMBINED LEGS W/ SYNCHRONIZED RECLINE & LEGS
31	X	X				X			TILT/RECLINE/ARTICULATING FOOT PLATFORM W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
32	X	X			X			X	TILT/RECLINE/COMBINED LEGS,/BACK W/ SYNCHRONIZED RECLINE & LEGS
33	X	X				X		X	TILT/RECLINE/ARTICULATING FOOT PLATFORM/BACK W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
34	X	X	X						TILT/ RECLINE/ & ELEVATING SEAT
35	X	X	X	X			X		TILT/ RECLINE/ ELEVATING SEAT/ INDEPENDENT LEGS W/ SYNCHRONIZED RECLINE & LEGS
36	X	X	X		X				TILT /RECLINE/ELEVATING SEAT/ COMBINED LEGS W/ SYNCHRONIZED RECLINE & LEGS
37	X	X	X			X			TILT/RECLINE/ELEVATING SEAT/ ARTICULATING FOOT PLATFORM W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
38	X	X	X		X			X	TILT /RECLINE/ELEVATING SEAT/ COMBINED LEGS/BACK W/ SYNCHRONIZED RECLINE & LEGS
39	X	X	X			X		X	TILT/RECLINE/ELEVATING SEAT/ ARTICULATING FOOT PLATFORM/ BACK W/ SYNCHRONIZED RECLINE & FOOT PLATFORM
LIFT & TILT CONFIGURATIONS (AAM)									
40	X		X						LIFT/ TILT
41	X		X	X			X		LIFT/TILT/INDEPENDENT LEGS
42	X		X		X				LIFT/TILT/COMBINED LEGS
43	X		X			X			LIFT/TILT/ARTICULATING FOOT PLATFORM
44	X	X	X						LIFT/TILT/RECLINE
45	X	X	X	X			X		LIFT/TILT/RECLINE/INDEPENDENT LEGS W/ SYNCHRONIZED RECLINE & LEGS
46	X	X	X		X				LIFT/TILT/RECLINE/COMBINED LEGS W/ SYNCHRONIZED RECLINE & LEGS
47	X	X	X			X			LIFT/TILT/RECLINE/ARTICULATING FOOT PLATFORM W/ SYNCHRONIZED RECLINE & FOOT PLATFORM