



Say hello to Firefly





GottaGo

Toileting Made Easy, Home and Away

An innovative, portable toilet seat which makes toilet training easier for families of children with special needs.



Product Features

- Suitable for children aged 2-9 years with a maximum user weight of 30kg/66lbs.
- GottaGo frame has 3 height settings: potty height, standard toilet height and disabled toilet height.
- Warm-touch, wipe clean materials guarantee a comfortable and snug fit.
- Easy to set up and put away.
- Specially tailored toilet training guide to support and motivate families throughout their toilet and potty training journey.



Technical Details		
	Size 1	Size 2
Age (approx.)	2-5 years	4-9 years
Hip Width	22.5cm	27cm
Minimum Seat Depth*	22cm	29cm
Backrest incline angle (forward)	8°	25°
Thigh angle (off horizontal)	40°	40°
Backrest height	24cm	36cm
Weight (seat and frame)	5.8kg	6.1kg

* For boys sizing it is recommended to choose the larger size providing they meet the minimum seat depth measurement.





GoTo Seat

One Little Seat, Countless Possibilities

A lightweight, portable and stylish postural support seat that enhances children's lives through participation.



Product Features

- Suitable for attachment to most standard chairs in the home, restaurants and classroom by using the adjustable attachment straps. Can also be attached to shopping trolleys.
- Adjustable headrest and laterals that adjust for growth.
- Standard or Advanced Headrest available.
- Pressure relief cushion as standard.
- Floorsitter has 26 recline positions, from upright (106°) to right back (140°).



Available in 2 fabric colours
and a vinyl version



Technical Details

	Size 1	Size 2
Age	1-3	3-8
Seat Base to Top of Headrest	380mm-540mm (14"-21")	510mm-680mm (20"-27")
Lateral Support Width	190mm-260mm (7.5"-10")	210mm-300mm (8"-12")
Seat Depth	175mm (7")	244mm (9.5")
Max. User Weight	15kg (33lbs)*	30kg (66lbs)**

*If a child is above 13kg/28.5lbs GoTo Seat cannot be used in a shopping trolley.

**Size 2 GoTo Seats are too large for use in shopping trolleys.

firefly Upsee

Little Steps to an Unexplored World

A mobility harness designed to allow children with special needs to experience upright movement like anyone else.



Product Features

- Facilitates gradual weight-bearing and gives children with motor impairment the closeness and confidence needed to stand and move with the help of an adult.
- Creates opportunity for social and emotional development as users can interact with peers at eye-level.
- Comprised of three parts - child harness, adult belt and double sandal.
- Durable, machine washable and breathable fabrics.
- Two colour options available for each size - please visit our website for further details.



Available in 4 fabric colours

Technical Details

	Extra Small	Small	Medium	Large
Age	1-2	2-4	4-6	6-8
Max. User Weight	15kg (33lbs)	15kg (33lbs)	20kg (44lbs)	25kg (55lbs)
Shoulders to Hip	280mm-320mm (11"-13.5")	310mm-360mm (12"-14")	350mm-390mm (14"-15.5")	380mm-420mm (15"-16.5")
Chest Circumference	500mm-570mm (19"-22")	520mm-600mm (20"-24")	570mm-650mm (22"-26")	630mm-710mm (25"-28")
Lower Hip Circumference	500mm-570mm (19"-22")	520mm-600mm (20"-24")	580mm-670mm (23"-26")	640mm-730mm (25"-29")



firefly Scooot

Stylish Wheels for Little Explorers

A 4-in-1 sit and ride which gives children the freedom to explore and have fun in style.



Product Features

- Scooot has four configurations: Scoot, Ride, Pooosh and Crawl. Different configurations for different abilities and different ways to play.
- For added support in Scooot, the Advanced Backrest Cushion helps. It's an optional extra that gives support to the child's trunk. The laterals move up and down to give support where it's needed. Or remove them and just have a comfy cushioned backrest. The Advanced Backrest Cushion is easy to attach, wash, and can be bought separately.
- One size caters for children aged 2-6- the footrest adjusts to suit the child's leg-length as they grow.



Technical Details	
Size	One size
Age	2-6
Max. User Weight	22kg (48.5lbs)
Seat Depth	400mm (15.5")
Seat Width	230mm (9")
Hip to end of footplate	520mm-630mm (20.5"-24.5")
Backrest Width	260mm (10.5")
Backrest Height	300mm (12")
Wheel Arch Width	230mm (9")
Pooosh Pole	800mm-890mm



Four Configurations

- **Crawl** allows kids to move and explore on their tummy, helping to strengthen their upper body and tummy muscles as they lift their heads and shoulders.
- **Scoot** lets kids sit upright to move around using their feet. Supported by a lap belt and backrest, the child's trunk and pelvis are encouraged to stabilise, freeing up their hands to lift and handle objects, play with family, and explore their surroundings.
- **Ride** lets kids propel themselves around the floor by moving the wheels with their hands. Having fun in this function can build upper body strength, whilst exploration and participation can improve spatial awareness as well as cognitive, sensory and communication abilities.
- **Pooosh** lets somebody else do the hard work. If the child isn't ready to propel themselves yet (or they're just tired), attach the Pooosh handle. It helps get little ones used to the Scooot. And it means parents and therapists can help them follow the fun at home, school or therapy sessions.



firefly Playpak

Fun Therapy for Little Learners

A portable, versatile early intervention therapy kit which provides children with the building blocks they need to help develop key skills.



Product Features

- Lightweight and portable.
- More than just practical- Playpak looks great with its colourful, fun, child-friendly design.
- Removable fabric can be washed at 40°C and tumble dried at a low temperature.
- Online Therapy Booklet, featuring expert advice and instructions to help develop physical, cognitive, sensory and communication skills.
- Lightweight bag containing an entire activity centre, including rolls, wedges and supports (designed by clinical experts) that stick together in dozens of different ways to practice the early positions which form the building blocks for later abilities.



Helps develop 5 key skills:

- Back Lying (Supine)
- Tummy Lying (Prone)
- Floor Sitting (Long Sitting)
- Side Lying
- Four Point Kneeling



firefly WeGo

Everyday Adventures in Style

A stylish, lightweight pushchair which is perfect for everyday use by active families and little ones who like to be on the move.



Product Features

- Compatible with the Firefly GoTo Seat.*
- Lightweight and easy to fold.
- Accessories included as standard- sun canopy, foot muff, rain cover and basket.
- Recline angle up to 140 degrees.
- Maximum user weight of 30kgs.*
- Wheels suitable for rough terrain.
- Small footprint so can fit in most car boots.

* Please check sizing guide for exact measurements



Technical Details	
Footprint Length (Collapsed)	116cm
Footprint Height (Collapsed)	38cm
WeGo Height	122cm
Seat Depth	26cm
Seat Height	50cm
Seat Width	35cm
Footrest Shortest Setting	25cm
Footrest Longest Setting	33cm
WeGo Weight	11kg
WeGo Weight with GoTo Seat	13kg

Please note - Not all children using the GoTo Seat Size 2 will be suitable for the WeGo. Please refer to our sizing guide.



Part of the
LECKEY
family circle

firefly **Splashy**

Water Adventures for Little Splashers

A lightweight, portable multipurpose bath seat that can be used for both bathing and other fun water play activities.



Product Features

- Splashy is designed to sit low in the bath to allow the child to benefit from being immersed in the water.
- Great aesthetics via bright colour options and a modern design.
- Splashy adapts to the child's postural support needs, with the option of a 5-point harness or the 3-point pelvic support. There are four detachable bumpers to add extra support wherever it's needed. The harness and foam supports attach easily to the seat base using toggles.
- Warm-touch foam feels soft whether it's wet or dry.
- Splashy's base unit lets you switch easily between 26 recline positions, from upright (106 degrees) to full recline (140 degrees).
- Lightweight, detachable parts make it easy to clean.



Available in 3 colours



Technical Details

Age	1-8
Backrest Height	610mm (24")
Backrest Width	385mm (15")
Backrest Depth	300mm (12")
Seat Width	350mm (14")
Floorsitter Width	420mm (16.5")
Floorsitter Length	630mm (24")
Max Weight	30kg (66lbs)



firefly **Splashy Big**

Water Adventures for Big Splashers

With all of the features and benefits of Splashy, the Splashy Big opens up many more opportunities for 8-16 year olds who want to participate in fun activities with friends and family.

Product Features

- Suitable for 8-16 years with a max user weight of 55KG/122lbs.
- Splashy Big is lightweight and super easy to maintain with its warm touch and quick drying covers. Cleaning it down after use is a breeze.
- Splashy Big's floorsitter has 26 recline options plus a new stabiliser feature for use outside of the bath.



Available in 2 colour combinations

Technical Details

Seat Backrest Height	720mm [28.5"]
Seat Backrest Width	450mm [17.5"]
Seat Base Depth	400mm [15.5"]
Seat Base Width	450mm [17.5"]
Floorsitter Width	420mm [16.5"]
Floorsitter Width with Stabiliser	665mm [26"]
Floorsitter Length	795mm [31"]
Maximum Weight	55kg [121lbs]
Age Range	8-16 years old



Notes