

PRODUCT SPECIFICATION

BEEASURE VIBE ASTM LEVEL 3 FACE MASKS

(ITEM #: BE2500/BE2510/BE2530)

1 PERFORMANCE

- Bacterial Filtration Efficiency(BFE)** : >99%
- Particle Filtration Efficiency(PEF)** : >99%
- Delta-P (ΔP)** : <5.0 mmH₂O/cm²
- ***In compliance to EN 14683(TYPE II) and ASTM Level 3 standards***

2 Finished Face Mask Size : 175*95(±5)mm

3 MATERIALS

- A. Observe Layer** : 100%PP Printed Spunbond (hydrophobic)
Basis weight: 28gsm
175mm±5(L)X200mm±5mm(W)
- B. Second Layer** : 100%PP Spunbond (hydrophobic)
Basis weight: 28gsm (white or carbon)
175mm±5(L)X178mm±5mm(W)
- C. Third Layer** : 100%PP Meltblown (hydrophobic)
Basis weight: 30gsm, white
175mm±5(L)X178mm±5mm(W)
- D. Inner Layer** : 100%PP Spunbond (hydrophobic)
Basis weight: 15gsm, white
175mm±5(L)X200mm±5mm(W)
- E. Ear loop** : Latex Free Elastic Loop
Spandex polyurethane coated with nylon thread
165mm±5(L)X3mm±0.5(W)
- F. Nose Bridge** : Full Aluminum nosepiece
Top:110mm±5mm(L)X3mm(W)X0.4mm(T)
Bottom: 90mm±5mm(L)X3mm(W)X0.4mm(T)
- G. Nose Flap** : Film Nonwoven
Basis weight: 40gsm
175mm±5 mm(L)X18mm±3mm & 42mm±3mm (W)

4 OBSERVE LAYER COLOUR : Luminous Blue, Dazzling Green, Glamorous Pink

5 PACKING : 50 pcs per dispenser box
8 boxes per carton

