

LECKEY BeMe

Chassis Size Guide

	Small	Medium Ind Gas	Medium Ind Hydraulic	Medium Ind/Out	Large Hydraulic	Large Power
Footprint Width [mm]	535	690	690	740	800	800
Footprint Length [mm]	635	840	840	880	770/850	770
Tilt Range [degrees]	- 10 to 30	- 10 to 25	- 10 to 25	- 10 to 25	- 10 to 30	- 10 to 30
Floor to Seat - Minimum	270	340	325	380	460	460
Floor to Seat - Maximum	650	655	615	670	640	640
Max Weight Limit [kg]**	35	50	50	50	70	70

**System max user weight should be lesser of seat, chassis and backrest