



## DELUXE MEDIUM-SIZED MOBILITY SCOOTER

Brand new forged frame and novel appearance design makes for the most popular mid-sized CTM model, with comfortable seat, pneumatic tires, rear suspension system and larger battery for enhanced stability, on top of LED full light control system for clear night lighting suitable for daily outdoor use.



# HS-520

# Features



Easy-To-Use Control Panel



Multi-Angle Adjustable Steering



Adjustable front and rear Seat



Revolving Tilting Seat



Adjustable Armrests



Charging Port At Tiller Head



Anti Tipping Safety Wheels

# Options



Cup Holder



Rear Basket



Walker Holder



Mobile phone holder

# Colors



Silver

Burgundy

Fine White



## CHIEN TI ENTERPRISE CO., LTD.

No. 13, Lane 227, Fu Ying Rd., Hsin Chuang District, New Taipei City, Taiwan  
 Tel : 886-2-2903-2987 Fax : 886-2-2904-1170 · 2903-8807  
 E-Mail:sales@chienti.com.tw <http://www.chienti.com.tw>

Overall Length	1200 mm / 47.2"
Overall Width	625 mm / 24.6"
Overall Height	1155 mm / 45.5"
Front Wheels	290 mm / 11"
Rear Wheels	290 mm / 11"
Weight W/ Batteries	94 kg / 208 lbs
Max. Speed	12 kmph / 7.5 mph
Weight Capacity	136 kg / 300 lbs
Ground Clearance	65 mm / 2.6"
Grade Climbable	10 degree
Curb Climbable	70 mm / 2.8"
Turning Radius	1400 mm / 55.1"
Suspension	Front & Rear
Brake	Electro-Mechanical
Seat Type	Swivel Mid-Back W/ Headrest & Seat Sliding Mechanism
Seat Width	465 mm / 18.3"
Motor Size	700W, 4500 r.p.m
Battery Size	(2) 12V. 50Ah
Weight of Battery Pack	30.5 kg / 67.2 lbs
Travel Range	39 km / 24.2 Miles
Battery Charger	5A Off Board
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Knob

\*Subject to change without notice.