



TOP CLASS COCKPIT LARGE-SIZED MOBILITY SCOOTER

CTM's largest, car-shaped, medical-grade mobility scooter with a professional waterproof cockpit, durable and clear windshield with wiper, detachable side doors to keep warm and comfortable all seasons, along with optional lighting and rearview camera for large body size and long distance needs.



HS-928C

Features



Front Windshield Wiper



Double Rear Mirrors



Full Light Control System



Standard Delta Handle Bar



LED Control Panel



Windshield Wiper & Temperature Panel



Air Vent Handle Bar

Options



LED Reading Light



Heating System



Rear View Camera

Colors



Silver

Overall Length	1600 mm / 63"
Overall Width	800 mm / 31.5" (Car body)
Overall Height	1620 mm / 63.8"
Front Wheels	380 mm / 15"
Rear Wheels	380 mm / 15"
Weight W/ Batteries	229 kg / 505 lbs
Max. Speed	10 kmph / 6.2 mph
Weight Capacity	200 kg / 441 lbs
Ground Clearance	110 mm / 4.3"
Grade Climbable	10 degree
Curb Climbable	100 mm / 3.9"
Turning Radius	1860 mm / 73.2"
Suspension	Front & Rear
Brake	Hand Brake & Electro-Mechanical
Seat Type	Swivel Captain Seat With Seat Sliding Mechanism
Seat Width	559 mm / 22"
Motor Size	1000W, 3200 r.p.m.
Battery Size	(2) 12V. 100Ah
Weight of Battery Pack	65.3 kg / 144 lbs
Travel Range	44 km / 27 Miles
Battery Charger	12A Off Board
Electronics	On / Off Key Switch, Battery Level Indicator, Speed Control Button

*Subject to change without notice.



CHIEN TI ENTERPRISE CO., LTD.

No. 13, Lane 227, Fu Ying Rd., Hsin Chuang District, New Taipei City, Taiwan
 Tel : 886-2-2903-2987 Fax : 886-2-2904-1170 · 2903-8807
 E-Mail:sales@chienti.com.tw http://www.chienti.com.tw